



„You can also integrate a 3 cm thick and 60 cm high skai-padding using the conversion kit available.“



Trendelenburg-funktion (Option)



Head & foot element manual or electric adjustment (Option)



Electric height adjustment



Thanks to the construction, a skai-padding also be subsequently adapted by a conversion kit.

Kathrin Klein
(Sales Department)



EN 60601-2-52:2010
DIN EN 32623:2010

Medical Devices Act
CE-Markering
Declaration of Conformity

LUKE

Why you should purchase Luke:

- Easy door handling with patented child proof lock
- Can be combined with all bed frames
- Optimum adjustment to correct care height possible, even for tall people
- The wooden surface is treated suitably for children (DIN 71-3 Toy Safety)
- Body weights of up to 150 kg possible
- Extensive accessories
- Professional construction style, high stability and reliability
- Available in various bed sizes (up to 180 x 220 cm)
- Quick delivery time
- Humidity protection
- Barred side height of 60 cm minimum (DIN/EN-Regulation) available!



Choice of colours
(see on page 71)



ID: 20633
4 winged door and back element with perspex

You can also integrate a 3 cm thick and 60 cm high skai-padding using the conversion kit available. The Skai padding is permanently mounted on a baseplate and can thus be integrated into the bed.

The bed is adapted for this so that all its functions continue to be available without restrictions.

technical data

lying area	external dimension	Article-no.	
90 x 200 cm	99 x 223 cm	92.32.009.L	
90 x 170 cm	99 x 193 cm	92.32.009.M	
doors			
folding door to open the complete front		99.703	•
folding door on both sides		99.750	+
locking mechanism			
patented lock IV		99.034A	•
lying height (top of slatted bed base)			
32-82 cm electric			
barred side height	overall height		
60 cm	99-149 cm		
slatted bed base			
no adjustment (fixed)			•
manual adjustment		99.030	+
electric adjustment		99.200	+
Trendelenburg-function (sloping position)		99.211	+
running gear			
four castor running gear with total locking mechanism			•
kind of wood			
beechwood			•
mattresses			
polyfoam mattress		99.031.4	•
interior spring mattress		99.031.2	+
reha-star mattress plus (up to stage 2)		99.031.3.PLUS	+
reha-star mattress II (up to stage 4)		99.031.3.ZWEI	+
incontinence cover		99.031.3ZA50	+

• = serial n.a.c. = no additional charge o.r. = on request + = Accessories

technical data

coloured features		
natural	90.100	•
yellow	90.128	n.a.c.
orange	90.204	n.a.c.
blue-purple	90.405	n.a.c.
dark-purple	90.406	n.a.c.
blue	90.515	n.a.c.
green	90.618	n.a.c.
equipment		
children's name on side of bed		n.a.c.
perspex instead of sticks		o.r.
external lock on side of bed	99.706	+
raising aid	99.802	+
mattress holder	99.891	+
protection for the bed base in case of kicks	99.892	+
wood care set	99.880	+
padding		
skai-padding	99.051LUKE90200 / 99.051LUKE90170	+
skai-padding with 4 vinyl windows	99.051LUKEVF90200 / 90170	+
retrofitkit: fixed skai-padding	99.051FESTLUKEUS90200 / 90170	+
fixed skai-padding (in order)	99.051FESTLUKE90200 / 90170	+
fixed skai-padding with 4 vinyl windows	99.051FESTLUKEVF90200/90170	+

• = serial n.a.c. = no additional charge o.r. = on request + = Accessories

LUKE

Fixed Skai-padding

More information about the different choices of padding please see page 61.

✓
upgradeable
fixed
skai-padding
(3 cm)



Locking of doors:



Patented lock IV

Integrated vinyl windows:

The integrated vinyl windows fit flush with the firm Skai padding and thus offer a similar padding effect to the permanently integrated foam padding. In addition it is possible to see in and out, which increases the scope for supervised use and the child's well-being.



LUKE

Remarkable features of Knut:

- Barred side height 60 cm
- Lying height 32 cm - 82 cm
- Trendelenburg-function (=Sloping position) (Option)
 - Perspex (Option)
- Apapted for fixed skai-padding (3 cm), retrofit kit is possible




ID: 19417



Protection for the bed base
in case of kicks

Our skai-padding is in
many colours available!
You can find a pre-selection on
our website at:
www.savi.de



Thanks to individual production, we can offer you the right cushioning for each specific need.

1. Bound skai-padding

Our skai-padding variants includes all listed below!

Version 1:

Skai-padding with Velcro-fixing, height adjustmend possible!
The Velcro webs are attached to the sticks.



Version 2:

The skai-padding is using cords on the side Sticks attached. Due to the positioning of the cords always flexible giben the optimal fixation.



Version 3:

Skai-padding for fixing with adhesive Velcro.



2. Motifpadding

Our motive padding is easy with cords at the bars of the bed fixed.



Greta



Paul



Paula

3. Fixed skai-padding

Fixed skaipadding, padding is screwed at the cot, cot-dimensions adapted so that you can use the adjustment. Attention: it's not possible to retrofit this version!



4. Fixed skai-padding with vinyl windows

The skaipadding is framing a vinyl-window. You and your child have a perfect view.



Our skai-padding is in many colours available!
You can find a preselection on our website at:
www.savi.de/Skaipadding



technical data

padding	
skai-padding 50 cm	99.051
skai-padding 100 cm	99.051.100
skai-padding with vinyl windows	
Motifpadding Greta, 50 cm	99.050GRETA
Motifpadding Paul, 50 cm	99.050PAUL
Motifpadding Paula, 50 cm	99.050PAULA
fixed skai-padding 50 cm	99.051.FEST
fixed skai-padding with vinyl windows	99.051.FEST
egde padding	99.052

thickness skai-padding 3 mm

5. Edge padding

The padding around the edges to protect against injury.



6. Fixed skai-padding with edge padding

Combine the fixed skai-padding with the edge padding.



7. Sticks- and Frame-skaipadding


Example pictures: please see page 33. All our children's care beds are available with special versions of skai-padding.

ACCESSORIES

Equipment, which makes nursing easier




Wood care set



Easy-close-mechanism




Touch protection




Cut out




Patented lock



Patented lock II




Patented lock III




Patented lock IV




Electric adjustable side rail



270° door



Removable corner posts




Protection for the bed base

ACCESSORIES



Children's name




Coloured sticks




Motor cover




Drawers



External lock on side of bed




Neet-roof



Feed tube hole




Feed tube hole



Entrance help



Raising aid



Play trainer wood



Play trainer metal

Create your own children's bed



Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

65 x 155 cm 70 x 140 cm 90 x 170 cm
90 x 200 cm 100 x 200 cm other size: _____

2. Which barred side do you need?

Height (cm): _____

3. Do you use a patient-lifter?

Yes No

4. Which lying height (top of slatted base) do you need?

40 cm 50 cm 60 cm 70 cm other height: _____

5. Do you need an electric height adjustment?

Yes No

6. Which slatted bed base do you need?

no adjustment manual adjustment protection for the bed base
electric adjustment Trendelenburg function

7. What kind of doors do you want?

2-winged door 4-winged door 4-winged sliding door 270° door
adjustable side rail electric adjustable side rail

8. How many doors do you need?

front side front and back sides head and footelement

9. Locking mechanism?

patented lock I patented lock II patented lock III patented lock IV

10. Do you need perspex?

Yes Where? _____ No

11. Kind of wood?

beech pine

12. Colour feature?

natural yellow orange blue-purple dark-purple
blue green white stained white lacquered other colour: _____

coloured sticks children's name: _____

13. Running gear?

wheels without wheels

14. Do you want containers?

Yes No

15. Do you want a mattress?

Yes No

polyfoam mattress interior spring mattress cold foam mattress

Reha-Star II Reha-Star Plus intercontinence cover

16. Do you need a padding?

Yes No

skai-padding fixed skai-padding skai-padding for the top edge of the cot

motif-padding: Paula Paul Greta

height of the padding: _____ colour of the padding: _____

17. Do you need Vinyl-windows?

Yes Where? _____ No

18. Equipment?

touch-protection net roof raising aid motor cover (Knut)

mattress-support side up aid removeable posts (Niklas)

- Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

- Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, data

stamp, sign



SAVI GmbH
Reha Solutions
Lise-Meitner-Str. 10
D-48691 Vreden

Postfach 1123
D-48684 Vreden
Fon: +49 (0) 25 64 - 93 09-70
Fax: +49 (0) 25 64 - 93 09-60

Very easy: Fax your order to
+49 25 64 - 93 09 - 60
or write an e-mail to
sales@savi.de
(Don't forget your contact details)

Or use our online-form at:
www.savi.de

www.savi.de
info@savi.de

Facebook:
SAVI.GmbH

SAVI APP



SAVI GmbH
Reha Solutions
Lise-Meitner-Straße 10
D-48691 Vreden

PO-Box 1123
D-48684 Vreden
Fon: +49(0)25 64 -93 09-70
Fax: +49(0)25 64 -93 09-60

please scan QR-Code



Your dealer: