



„For me as a purchaser it’s especially the priceperformance ratio that counts.  
And SAVI offers a great ratio.”

Thomas Mikosch  
(Kid(z) & Co.)



Electric height adjustment



Head & foot element manual or electric adjustment (Option)



Niklas is designed with extreme flexibility and functionality of the barred sides in mind.

The bed is accessible from all sides because the sides can be lowered independently of each other.

Top safety paired with unrestricted freedom – you can choose.



(Option)

EN 60601-2-52:2010  
DIN EN 32623:2010

Medical Devices Act  
CE-Markering  
Declaration of Conformity

# NIKLAS




Perspex + sticks ⓘ  
+ colour features deep-purple



All perspex  
+ electric adjustable sides  
+ colour features deep-purple




All perspex  
+ electric adjustable sides  
+ colour features deep-purple



Choice of colours  
(see on page 71)

## Why you should purchase Niklas:

- Easy door handling with patented child proof lock
- Can be combined with all bed frames
- Optimum adjustment to correct care height possible, even for tall people
- The wooden surface is treated suitably for children (DIN 71-3 Toy Safety)
- Body weights of up to 150 kg possible
- Extensive accessories
- Professional construction style, high stability and reliability
- Available in various bed sizes (up to 180 x 220 cm)
- Quick delivery time
- Humidity protection
- Barred side height of 60 cm minimum (DIN/EN-Regulation) available!



Three kinds of our motifpadding are available. You will find more informations on our website: [www.savi.de](http://www.savi.de)

## \* Indication

The electric height adjustment 47-87 cm is only in combination with a barred side height of 40 cm + 20 cm high slip-possible.

## technical data

| lying area   | external dimension | Article-no.           |                 |
|--|--------------------|-----------------------|-----------------|
| 100 x 200 cm   | 115 x 215 cm       | 131.50.006            |                 |
| 90 x 200 cm  | 105 x 215 cm       | 131.50.006.L          |                 |
| 90 x 170 cm  | 105 x 185 cm       | 131.50.006.M          |                 |
| 65 x 155 cm  | 80 x 170 cm        | 131.50.006.S          |                 |
| <b>doors</b>   |                    |                       |                 |
| all sides adjustable in steps (lowerable)  |                    |                       | •               |
| <b>lying height (top of slatted bed base)</b>  |                    |                       |                 |
| 76-116 cm electric height adjustment   |                    |                       | •               |
| 47- 87 cm electric height adjustment (only in combination with 99.737 possible, please see on page 18) |                    |                       | o.r..           |
| <b>barred side height</b>  |                    | <b>overall height</b> |                 |
| 60 cm  | 140 cm             |                       | •               |
| 40 cm + 20 cm *  | 155 cm             | 99.737                | +               |
| <b>slatted bed base</b>  |                    |                       |                 |
| no adjustment (fixed)  |                    |                       | +               |
| manual adjustment  |                    |                       | 99.030 •        |
| electric adjustment  |                    |                       | 99.201 +        |
| Trendelenburg-function A (sloping position)  |                    |                       | 99.212 +        |
| Trendelenburg-function B (sloping position)  |                    |                       | 99.213 +        |
| <b>running gear</b>  |                    |                       |                 |
| four castor running gear with total locking mechanism  |                    |                       | •               |
| designer wooden castors  |                    |                       | 99.885 n.a.c.   |
| <b>kind of wood</b>  |                    |                       |                 |
| beechwood  |                    |                       | •               |
| <b>mattresses</b>  |                    |                       |                 |
| polyfoam mattress  |                    |                       | •               |
| interior spring mattress   |                    |                       | 99.031.2 +      |
| reha-star mattress plus (up to stage 2)  |                    |                       | 99.031.3.PLUS + |
| reha-star mattress II (up to stage 4)  |                    |                       | 99.031.3.ZWEI + |
| incontinence cover   |                    |                       | 99.031.3ZA50 +  |

• = serial    n.a.c. = no additional charge    o.r. = on request    + = Accessories

## technical data

| coloured features                            | Article-no.    |        |
|--|----------------|--------|
| natural                                      | 90.100         | n.a.c. |
| yellow                                       | 90.128         | n.a.c. |
| orange                                       | 90.204         | n.a.c. |
| blue-purple                                  | 90.405         | n.a.c. |
| dark-purple                                  | 90.406         | •      |
| blue   | 90.515         | n.a.c. |
| green  | 90.618         | n.a.c. |
| <b>equipment</b>                             |                |        |
| children's name on side of bed               |                | n.a.c. |
| coloured sticks                              |                | n.a.c. |
| perspex + sticks (as shown) ⓘ                | 99.810LEXANNIK | +      |
| 4 removable corner posts                     | 99.889         | +      |
| only 1 side adjustable, 3 sides fixed        | 99.730         | +      |
| only 2 sides adjustable, 2 sides fixed       | 99.731         | +      |
| 1 electric adjustable side rail              | 99.735         | +      |
| 2 electric adjustable side rails             | 99.736         | +      |
| feed tube hole                               | 99.840         | o.r..  |
| raising aid                                  | 99.802         | +      |
| sit up aid                                   | 99.820         | +      |
| adaption of the bed for sloping ceiling      |                | o.r.   |
| 2 drawers                                    | 99.040         | +      |
| mattress holder                              | 99.891         | +      |
| protection for the bed base in case of kicks | 99.892         | +      |
| wood care set                                | 99.880         | +      |
| <b>padding</b>                               |                |        |
| motifpadding Greta, 50 cm                    | 99.050GRETA    | +      |
| motifpadding Paul, 50 cm                     | 99.050PAUL     | +      |
| motifpadding Paula, 50 cm                    | 99.050PAULA    | +      |
| skai-padding 50 cm height                    | 99.051NIKLAS   | +      |
| fixed skai-padding                           | 99.051.FEST    | o.r.   |
| fixed skai-padding with vinyl windows        | 99.051.FEST    | o.r.   |
| egde padding                                 | 99.052         | +      |

• = serial    n.a.c. = no additional charge    o.r. = on request    + = Accessories

# NIKLAS

## Remarkable features of Niklas:

- Barred side height 60 cm, all side-rails adjustable
- Trendelenburgfunktion (Option)
- Lying height 76 cm - 116 cm or 47 - 87 cm (low entrance)
- Electric adjustable sides (Option)
- Removeable corner posts (Option)
- Perspexs (as shown) or with sticks
- Choice of colours (see on page 71)



Perspexs + sticks  
+ low entrance 47 - 87 cm  
+ with top element



## Benefits of electric adjustable sides:

- Easy to use at the push of a button
- Lockable control element prevents unsupervised opening
- Can also be operated by people with disabilities
- Minimal operating strength required (pushbutton)
- Infinitely adjustable

## Recessed safety grip:



Recessed safety grip:  
(Basic)

- Easy to use manually
- The safety gap on both recesses prevents unsupervised opening
- Minimal operating strength required (slider buttons)
- Grid adjustment

## Removable corner posts:



# NIKLAS

## Fixed Skai-padding

More information about the different choices of padding please see page 61.



ID: 20841



ID: 16581



ID: 15368



ID: 17167



ID: 10352



ID: 15124



ID: 20631

# NIKLAS



ID: 10400  
Coloured sticks



ID: 16069  
Drawers

Trendelenburg-function A:



Inclination of the  
complete bed

Trendelenburg-function B:



Inclination of the  
slatted bed base



ID: 10682  
Canopy bed,  
total height 215 cm, entrance height 30cm



ID: 10394  
Adaption of the bed  
for sloping ceiling



ID: 20152  
Play trainer wood



ID: 19714  
Head- and footelement  
with perspex



ID: 19204  
Small castors



ID: 18692  
White stained



ID: 20214  
Closeable feed tube hole with sliding mechanism



ID: 18499  
Without castors



ID: 19923  
Drawers



Thanks to individual production, we can offer you the right cushioning for each specific need.

## 1. Bound skai-padding

Our skai-padding variants includes all listed below!

### Version 1:

Skai-padding with Velcro-fixing, height adjustment possible! The Velcro webs are attached to the sticks.



### Version 2:

The skai-padding is using cords on the side Sticks attached. Due to the positioning of the cords always flexible given the optimal fixation.



### Version 3:

Skai-padding for fixing with adhesive Velcro.



## 2. Motifpadding

Our motive padding is easy with cords at the bars of the bed fixed.



Greta



Paul



Paula

## 3. Fixed skai-padding

Fixed skai-padding, padding is screwed at the cot, cot-dimensions adapted so that you can use the adjustment. Attention: it's not possible to retrofit this version!




## 4. Fixed skai-padding with vinyl windows

The skai-padding is framing a vinyl-window. You and your child have a perfect view.



Our skai-padding is in many colours available! You can find a preselection on our website at: [www.savi.de/Skaipadding](http://www.savi.de/Skaipadding)



## technical data

|                                       |             |
|---------------------------------------|-------------|
| padding                               |             |
| skai-padding 50 cm                    | 99.051      |
| skai-padding 100 cm                   | 99.051.100  |
| skai-padding with vinyl windows       |             |
| Motifpadding Greta, 50 cm             | 99.050GRETA |
| Motifpadding Paul, 50 cm              | 99.050PAUL  |
| Motifpadding Paula, 50 cm             | 99.050PAULA |
| fixed skai-padding 50 cm              | 99.051.FEST |
| fixed skai-padding with vinyl windows | 99.051.FEST |
| egde padding                          | 99.052      |

thickness skai-padding 3 mm

## 5. Edge padding

The padding around the edges to protect against injury.



## 6. Fixed skai-padding with edge padding

Combine the fixed skai-padding with the edge padding.



## 7. Sticks- and Frame-skaipadding


Example pictures: please see page 33. All our children's care beds are available with special versions of skai-padding.

# ACCESSORIES

Equipment, which makes nursing easier




Wood care set



Easy-close-mechanism



Touch protection



Cut out



Patented lock



Patented lock II




Patented lock III




Patented lock IV




Electric adjustable side rail



270° door



Removable corner posts



Protection for the bed base

# ACCESSORIES



Children's name



Coloured sticks




Motor cover




Drawers



External lock on side of bed




Neet-roof



Feed tube hole



Feed tube hole



Entrance help



Raising aid



Play trainer wood



Play trainer metal

# Create your own children's bed



## Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

### 1. Lying area of the bed?

65 x 155 cm  70 x 140 cm  90 x 170 cm   
90 x 200 cm  100 x 200 cm  other size: \_\_\_\_\_

2. Which barred side do you need? Height (cm): \_\_\_\_\_

3. Do you use a patient-lifter? Yes  No

### 4. Which lying height (top of slatted base) do you need?

40 cm  50 cm  60 cm  70 cm  other height: \_\_\_\_\_

5. Do you need an electric height adjustment? Yes  No

### 6. Which slatted bed base do you need?

no adjustment  manual adjustment  protection for the bed base   
electric adjustment  Trendelenburg function

### 7. What kind of doors do you want?

2-winged door  4-winged door  4-winged sliding door  270° door   
adjustable side rail  electric adjustable side rail

8. How many doors do you need? front side  front and back sides  head and footelement

### 9. Locking mechanism?

patented lock I  patented lock II  patented lock III  patented lock IV

### 10. Do you need perspex?

Yes  Where? \_\_\_\_\_ No

### 11. Kind of wood?

beech  pine

### 12. Colour feature?

natural  yellow  orange  blue-purple  dark-purple   
blue  green  white stained  white lacquered  other colour: \_\_\_\_\_

coloured sticks  children's name: \_\_\_\_\_

### 13. Running gear?

wheels  without wheels

### 14. Do you want containers?

Yes  No

### 15. Do you want a mattress?

Yes  No

polyfoam mattress  interior spring mattress  cold foam mattress

Reha-Star II  Reha-Star Plus  intercontinence cover

### 16. Do you need a padding?

Yes  No

skai-padding  fixed skai-padding  skai-padding for the top edge of the cot

motif-padding: Paula  Paul  Greta

height of the padding: \_\_\_\_\_ colour of the padding: \_\_\_\_\_

### 17. Do you need Vinyl-windows?

Yes  Where? \_\_\_\_\_ No

### 18. Equipment?

touch-protection  net roof  raising aid  motor cover (Knut)

mattress-support  side up aid  removeable posts (Niklas)

- Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

- Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, data

stamp, sign



SAVI GmbH  
Reha Solutions  
Lise-Meitner-Str. 10  
D-48691 Vreden

Postfach 1123  
D-48684 Vreden  
Fon: +49 (0) 25 64 - 93 09-70  
Fax: +49 (0) 25 64 - 93 09-60

Very easy: Fax your order to  
+49 25 64 - 93 09 - 60  
or write an e-mail to  
sales@savi.de  
(Don't forget your contact details)

Or use our online-form at:  
**www.savi.de**

[www.savi.de](http://www.savi.de)  
[info@savi.de](mailto:info@savi.de)

Facebook:  
SAVI.GmbH

SAVI APP



SAVI GmbH  
Reha Solutions  
Lise-Meitner-Straße 10  
D-48691 Vreden

PO-Box 1123  
D-48684 Vreden  
Fon: +49(0)25 64 -93 09-70  
Fax: +49(0)25 64 -93 09-60

please scan QR-Code



Your dealer: