



„For us it was important, that SAVI sets patterns for ergonomics and ecology.“



Electric height adjustment



Head & foot element manual or electric adjustment (Option)



This model combines safety with comfort.

The electric lying-height adjustment is especially useful when it comes to the daily care of handicapped people.

Familie Boland
(Westfalen)



EN 60601-2-52:2010
DIN EN 32623:2010

Medical Devices Act
CE-Markering
Declaration of Conformity

SAVOIR-VIVRE



Barred side height 60 cm
+ colour features green



Barred side height 100 cm
+ colour features blue



Barred side height 125 cm
+ colour features orange

Why you should purchase Savoir-Vivre:

- Easy door handling with patented child proof lock
- Can be combined with all bed frames
- Optimum adjustment to correct care height possible, even for tall people
- The wooden surface is treated suitably for children (DIN 71-3 Toy Safety)
- Body weights of up to 150 kg possible
- Extensive accessories
- Professional construction style, high stability and reliability
- Available in various bed sizes (up to 180 x 220 cm)
- Quick delivery time
- Humidity protection
- Barred side height of 60 cm minimum (DIN/EN-Regulation) available!



Barred side height 150 cm
+ colour features deep-purple




Barred side height 170 cm
+ colour features blue-purple

technical data

| lying area | external dimension | Article-no. |
|---|-----------------------|---------------------|
| 100 x 200 cm | 113 x 206 cm | please specify size |
| 90 x 200 cm (Size L) | 103 x 206 cm | please specify size |
| 90 x 170 cm (Size M) | 103 x 176 cm | please specify size |
| 65 x 155 cm (Size S) | 78 x 161 cm | please specify size |
| doors | | |
| 4-winged door to open the complete front | | • |
| four-winged door on both sides | 99.750 | + |
| 270° door, door outwardly opening | 99.704 | + |
| 270° door, door outwardly opening on both sides | 99.704BEID | |
| locking mechanism | | |
| patented child-prooflock | | • |
| patented child-prooflock metal | 40100003 | + |
| patented lock II | 99.033 | + |
| patented lock III | 99.034 | + |
| patented lock IV | 99.034A | + |
| lying height (top of slatted bed base) | | |
| 40-80 cm electric height adjustment | | • |
| Barred side height | | |
| | Overall height | |
| 60/70 cm | 110 cm | 110.50.504.99101 |
| 100/110 cm | 150 cm | 150.50.504.99101 |
| | for size 90 x 170 cm | |
| 125/135 cm | 175 cm | 175.50.504.99101 |
| | for size 90 x 170 cm | |
| 150/160 cm | 200 cm | 200.50.504.99101 |
| 170/180 cm | 220 cm | 220.50.504.99101 |
| special sizes on inquire | | |
| slatted bed base | | |
| no adjustment (fixed) | | + |
| manual adjustment | 99.030 | • |
| electric adjustment | 99.200 | + |
| Trendelenburg-function B (sloping position) | 99.213 | + |
| running gear | | |
| four castor running gear with total locking mechanism | 99.887 | + |
| designer wooden castors | 99.886 | + |
| cut out on the front side (notify height please) | 99.830 | n.a.c. |
| kind of wood | | |
| pinewood | | n.a.c. |
| beechwood | | • |
| Mattresses | | |
| polyfoam mattress | | • |
| interior spring mattress | 99.031.2 | + |
| reha-star mattress plus (up to stage 2) | 99.031.3.PLUS | + |
| reha-star mattress II (up to stage 4) | 99.031.3.ZWEI | + |
| incontinence cover | 99.031.3ZA50 | + |

• = serial n.a.c. = no additional charge o.r. = on request + = Accessories



Our skai-padding is in many colours available!
You can find a pre-selection on our website at:
www.savi.de

technical data

| coloured features | | |
|--|-------------|--------|
| natural | 90.100 | n.a.c. |
| yellow | 90.128 | n.a.c. |
| orange | 90.204 | n.a.c. |
| blue-purple | 90.405 | n.a.c. |
| dark-purple | 90.406 | n.a.c. |
| blue | 90.515 | n.a.c. |
| green | 90.618 | • |
| equipment | | |
| children's name on side of bed | | n.a.c. |
| coloured sticks | | n.a.c. |
| perspex instead of sticks | | o.r. |
| feed tube hole | 99.840 | o.r. |
| external lock on side of bed | 99.706 | + |
| raising aid | 99.802 | + |
| sit up aid | 99.820 | + |
| neet-roof | 99.800 | + |
| adaption of the bed for sloping ceiling | | o.r. |
| mattress holder | 99.891 | + |
| protection for the bed base in case of kicks | 99.892 | + |
| wood care set | 99.880 | + |
| padding | | |
| Motifpadding Greta, 50 cm | 99.050GRETA | + |
| Motifpadding Paul, 50 cm | 99.050PAUL | + |
| Motifpadding Paula, 50 cm | 99.050PAULA | + |
| skai-padding 50 cm height | 99.051 | + |
| skai-padding 100 cm height | 99.051.100 | + |
| fixed skai-padding | 99.051.FEST | o.r. |
| fixed skai-padding with vinyl windows | 99.051.FEST | o.r. |
| egde padding | 99.052 | + |

• = serial n.a.c. = no additional charge o.r. = on request + = Accessories

Remarkable features of Savoir-Vivre:

- Barred side height 60 cm - 170 cm or special sizes on inquire
- Electric lying height 40 cm - 80 cm
- Four-winged sliding door
- Four-casters running gear (Option)
- Patented child-prooflock
- Perspex (Option)

Barred side height Door height Total height

- | | | |
|----------|----------|----------|
| • 60 cm | • 70 cm | • 110 cm |
| • 100 cm | • 110 cm | • 150 cm |
| • 125 cm | • 135 cm | • 175 cm |
| • 150 cm | • 160 cm | • 200 cm |
| • 170 cm | • 180 cm | • 220 cm |



Total height:

- ↑ 110 cm ↑ 150 cm ↑ 175 cm ↑ 200 cm ↑ 220 cm

Locking of doors:



External lock on side of bed
(Option)

Cut out:



Height: 8 cm (Option)
(special sizes on require)



ID: 18384
2 x patented lock IV



ID: 16662
All perspex



ID: 19735
4-caster running gear
with total lockable mechanism



ID: 18579
Stained white



ID: 17448
Coloured sticks



ID: 17488
Cut out



ID: 19059
Back element with wooden panels



ID: 18454
Drawers



ID: 19064
Head and foot element with wooden panels

SAVOIR-VIVRE

Fixed Skai-padding

More information about the different choices of padding please see page 61.



ID: 17150



ID: 18629



ID: 19833



ID: 17098



ID: 19140



ID: 19037



ID: 17082


fixed Skai-padding:



ID: 19778

Integrated vinyl windows:

The integrated vinyl windows fit flush with the firm Skai padding and thus offer a similar padding effect to the permanently integrated foam padding. In addition it is possible to see in and out, which increases the scope for supervised use and the child's well-being.



Our skai-padding is in many colours available! You can find a pre-selection on our website at: www.savi.de



ID: 19931



ID: 19931



ID: 18542



ID: 19593



ID: 17934

Thanks to individual production, we can offer you the right cushioning for each specific need.

1. Bound skai-padding

Our skai-padding variants includes all listed below!

Version 1:

Skai-padding with Velcro-fixing, height adjustment possible! The Velcro webs are attached to the sticks.



Version 2:

The skai-padding is using cords on the side Sticks attached. Due to the positioning of the cords always flexible given the optimal fixation.



Version 3:

Skai-padding for fixing with adhesive Velcro.



2. Motifpadding

Our motive padding is easy with cords at the bars of the bed fixed.



Greta



Paul



Paula

3. Fixed skai-padding

Fixed skai-padding, padding is screwed at the cot, cot-dimensions adapted so that you can use the adjustment. Attention: it's not possible to retrofit this version!




4. Fixed skai-padding with vinyl windows

The skai-padding is framing a vinyl-window. You and your child have a perfect view.



Our skai-padding is in many colours available! You can find a preselection on our website at: www.savi.de/Skaipadding



technical data

| | |
|---------------------------------------|-------------|
| padding | |
| skai-padding 50 cm | 99.051 |
| skai-padding 100 cm | 99.051.100 |
| skai-padding with vinyl windows | |
| Motifpadding Greta, 50 cm | 99.050GRETA |
| Motifpadding Paul, 50 cm | 99.050PAUL |
| Motifpadding Paula, 50 cm | 99.050PAULA |
| fixed skai-padding 50 cm | 99.051.FEST |
| fixed skai-padding with vinyl windows | 99.051.FEST |
| egde padding | 99.052 |

thickness skai-padding 3 mm

5. Edge padding

The padding around the edges to protect against injury.



6. Fixed skai-padding with edge padding

Combine the fixed skai-padding with the edge padding.



7. Sticks- and Frame-skaipadding


Example pictures: please see page 33. All our children's care beds are available with special versions of skai-padding.

ACCESSORIES

Equipment, which makes nursing easier




Wood care set



Easy-close-mechanism



Touch protection



Cut out



Patented lock



Patented lock II




Patented lock III




Patented lock IV




Electric adjustable side rail



270° door



Removable corner posts



Protection for the bed base

ACCESSORIES



Children's name



Coloured sticks




Motor cover




Drawers



External lock on side of bed




Neet-roof



Feed tube hole



Feed tube hole



Entrance help



Raising aid



Play trainer wood



Play trainer metal

Create your own children`s bed



Create your children`s bed for rehabilitation!

This is a list which you can use to create your children`s bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed`s.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

65 x 155 cm 70 x 140 cm 90 x 170 cm
90 x 200 cm 100 x 200 cm other size: _____

2. Which barred side do you need? Height (cm): _____

3. Do you use a patient-lifter? Yes No

4. Which lying height (top of slatted base) do you need?

40 cm 50 cm 60 cm 70 cm other height: _____

5. Do you need an electric height adjustment? Yes No

6. Which slatted bed base do you need?

no adjustment manual adjustment protection for the bed base
electric adjustment Trendelenburg function

7. What kind of doors do you want?

2-winged door 4-winged door 4-winged sliding door 270° door
adjustable side rail electric adjustable side rail

8. How many doors do you need? front side front and back sides head and footelement

9. Locking mechanism?

patented lock I patented lock II patented lock III patented lock IV

10. Do you need perspex?

Yes Where? _____ No

11. Kind of wood?

beech pine

12. Colour feature?

natural yellow orange blue-purple dark-purple
blue green white stained white lacquered other colour: _____

coloured sticks children`s name: _____

13. Running gear?

wheels without wheels

14. Do you want containers?

Yes No

15. Do you want a mattress?

Yes No

polyfoam mattress interior spring mattress cold foam mattress

Reha-Star II Reha-Star Plus intercontinence cover

16. Do you need a padding?

Yes No

skai-padding fixed skai-padding skai-padding for the top edge of the cot

motif-padding: Paula Paul Greta

height of the padding: _____ colour of the padding: _____

17. Do you need Vinyl-windows?

Yes Where? _____ No

18. Equipment?

touch-protection net roof raising aid motor cover (Knut)

mattress-support side up aid removeable posts (Niklas)

- Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

- Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, data

stamp, sign



SAVI GmbH
Reha Solutions
Lise-Meitner-Str. 10
D-48691 Vreden

Postfach 1123
D-48684 Vreden
Fon: +49 (0) 25 64 - 93 09-70
Fax: +49 (0) 25 64 - 93 09-60

Very easy: Fax your order to
+49 25 64 - 93 09 - 60
or write an e-mail to
sales@savi.de
(Don't forget your contact details)

Or use our online-form at:
www.savi.de

www.savi.de
info@savi.de

Facebook:
SAVI.GmbH

SAVI APP



SAVI GmbH
Reha Solutions
Lise-Meitner-Straße 10
D-48691 Vreden

PO-Box 1123
D-48684 Vreden
Fon: +49(0)25 64 -93 09-70
Fax: +49(0)25 64 -93 09-60

please scan QR-Code



Your dealer: