



„I've never had two customers with exactly the same requirements.  
SAVI could always provide an individual solution.“



Electric height  
adjustment



Head & foot  
element manual  
or electric  
adjustment  
(Option)



Thanks to the four-winged folding sliding door, the barred sides and the flexible height adjustment possibilities, this bed provides very good therapeutic support.

Due to the low entry height of a little about 30 cm even smaller children can get into the bed on their own.

Günter Bußmann

(Manager)



EN 60601-2-52:2010  
DIN EN 32623:2010

Medical Devices Act  
CE-Markering  
Declaration of Conformity

# T O M



Barred side height 113 cm  
+ low entrance height of 20 cm  
+ colour features blue-purple



Barred side height 135 cm  
+ wooden castors  
+ colour features blue-purple



Barred side height 72 cm  
+ wooden castors  
+ colour features blue-purple



Choice of colours  
(see on page 71)


## Why you should purchase Tom:

- Easy door handling with patented child proof lock
- Can be combined with all bed frames
- Optimum adjustment to correct care height possible, even for tall people
- The wooden surface is treated suitably for children (DIN 71-3 Toy Safety)
- Body weights of up to 150 kg possible
- Extensive accessories
- Professional construction style, high stability and reliability
- Available in various bed sizes (up to 180 x 220 cm)
- Quick delivery time
- Humidity protection
- Barred side height of 60 cm minimum (DIN/EN-Regulation) available!

## technical data

lying area	external dimension	Article-no.	
100 x 200 cm	115,5 x 223 cm	155.34.008	
90 x 200 cm	105,5 x 223 cm	155.34.008.L	
90 x 170 cm	105,5 x 193 cm	155.34.008.M	
65 x 155 cm	80,5 x 178 cm	155.34.008.S	
70 x 140 cm	85,5 x 163 cm	155.34.008.XS	
<b>doors</b>			
four-winged door to open the complete front			•
four-winged door on both sides		99.750	+
270° door, door outwardly opening		99.704	+
270° door, door outwardly opening on both sides		99.704BEID	+
<b>locking mechanism</b>			
patented child-prooflock			•
patented child-prooflock metal		40100003	+
patented lock II		99.033	+
patented lock III		99.034	+
patented lock IV		99.034A	+
<b>lying height (top of slatted bed base)</b>			
34-104 cm electric height adjustment			•
<b>barred side height</b>	<b>overall height</b>		
72 cm	114 cm		+
113 cm	155 cm		•
135 cm	177 cm		+
special sizes on inquire			+
<b>slatted bed base</b>			
no adjustment (fixed)			+
manual adjustment			•
electric adjustment		99.200	+
Trendelenburg-function (sloping position)		99.211	+
<b>running gear</b>			
four castor running gear with total locking mechanism			•
designer wooden castors		99.885	n.a.c.
without castors and low entrance height of 20 cm		99.883	+
small castors an low entrance height of 27 cm		99.881	+
<b>kind of wood</b>			
beechwood			•
<b>mattresses</b>			
polyfoam mattress			•
interior spring mattress		99.031.2	o.r.
reha-star mattress plus (up to stage 2)		99.031.3.PLUS	+
reha-star mattress II (up to stage 4)		99.031.3.ZWEI	+
incontinence cover		99.031.3ZA50	+

• = serial    n.a.c. = no additional charge    o.r. = on request    + = Accessories



Our skai-padding is in many colours available!  
You can find a pre-selection on our website at:  
[www.savi.de](http://www.savi.de)

## technical data

<b>coloured features</b>		
natural	90.100	n.a.c.
yellow	90.128	n.a.c.
orange	90.204	n.a.c.
blue-purple	90.405	•
dark-purple	90.406	n.a.c.
blue	90.515	n.a.c.
green	90.618	n.a.c.
<b>equipment</b>		
children's name on side of bed		n.a.c.
coloured sticks		n.a.c.
perspex instead of sticks		•
sticks instead of perspex	99.815	n.a.c.
touch protection for height adjustment	99.900	+
external lock on side of bed	99.706	+
feed tube hole	99.840	o.r.
raising aid	99.802	+
neet-roof	99.800	+
adaption of the bed for sloping ceiling		o.r.
mattress holder	99.891	+
protection for the bed base in case of kicks	99.892	+
wood care set	99.880	+
<b>padding</b>		
skai-padding 50 cm height	99.051	+
skai-padding 100 cm height	99.051.100	+
fixed skai-padding	99.051.FEST	o.r.
fixed skai-padding with 4 vinyl windows	99.051.FEST	o.r.
egde padding	99.052	+

• = serial    n.a.c. = no additional charge    o.r. = on request    + = Accessories

# T O M


## Four-winged sliding door :

Thanks to its folding-sliding mechanism, the door offers the maximum stability and precise operation in every position. The entire front side can be opened in this way. The design is a harmonious combination of quality and function. Naturally by request we can also offer customised additional doors, and custom-made sloping roofs are also not a problem for us.



ID: 18965

## Locking of doors:



Patented child-prooflock (Base)



Patented child-prooflock metal



Patented lock II



Patented lock III:  
The labyrinth guide prevents the child opening the bed by accident. This patented lock is keyless and has no magnets. Can also be retrofitted.



Patented lock IV:  
Here too the bed is opened by means of a labyrinth guide, but in addition this lock also operates against resistance. This patented lock is keyless, without magnets and can be retrofitted.



External lock on side of bed  
(increases the stability and provides more secure)

## Benefits of the 270° door:

- Easy to use manually
- Doors on the head and foot elements fold back completely
- 270° version can also be supplied with firm Skai-padding



ID: 20109



ID: 20009



ID: 15938

## Remarkable features of Tom:

- Barred side height 72 cm, barred side height 113 cm, barred side height 135 cm or special sizes on inquire
- Slatted bed base 34 cm - 104 cm or slatted bed base 20 cm - 90 cm (without wheels)
- Trendelenburg-function (Sloping position) (Option)
- Choose between perspex or sticks



ID: 18463

## Fixed Skai-padding

More information about the different choices of padding please see page 61.



ID: 18677



ID: 20048



ID: 13876



ID: 18574



ID: 18471

## Vinyl windows:

The integrated vinyl windows fit flush with the firm skai-padding and thus offer a similar padding effect to the permanently integrated foam padding. In addition it is possible to see in and out, which increases the scope for supervised use and the child's well-being.



ID: 18866



ID: 20018



ID: 20723



ID: 20324



ID: 14364  
White painted



ID: 13175  
Low entrance height  
without castors



ID: 20011  
Coloured sticks



ID: 18429  
Adaption for sloping ceiling



ID: 16472  
Children's name  
on side of bed



ID: 14280  
Barred side height  
72 cm



ID: 17767  
All sticks



ID: 14849  
Touch protection  
for height adjustment

Thanks to individual production, we can offer you the right cushioning for each specific need.

## 1. Bound skai-padding

Our skai-padding variants includes all listed below!

### Version 1:

Skai-padding with Velcro-fixing, height adjustmend possible!  
The Velcro webs are attached to the sticks.



### Version 2:

The skai-padding is using cords on the side Sticks attached. Due to the positioning of the cords always flexible giben the optimal fixation.



### Version 3:

Skai-padding for fixing with adhesive Velcro.



## 2. Motifpadding

Our motive padding is easy with cords at the bars of the bed fixed.



Greta



Paul



Paula

## 3. Fixed skai-padding

Fixed skaipadding, padding is screwed at the cot, cot-dimensions adapted so that you can use the adjustment. Attention: it's not possible to retrofit this version!



## 4. Fixed skai-padding with vinyl windows

The skaipadding is framing a vinyl-window. You and your child have a perfect view.



Our skai-padding is in many colours available!  
You can find a preselection on our website at:  
[www.savi.de/Skaipadding](http://www.savi.de/Skaipadding)

## technical data

padding	
skai-padding 50 cm	99.051
skai-padding 100 cm	99.051.100
skai-padding with vinyl windows	
Motifpadding Greta, 50 cm	99.050GRETA
Motifpadding Paul, 50 cm	99.050PAUL
Motifpadding Paula, 50 cm	99.050PAULA
fixed skai-padding 50 cm	99.051.FEST
fixed skai-padding with vinyl windows	99.051.FEST
egde padding	99.052

thickness skai-padding 3 mm

## 5. Edge padding

The padding around the edges to protect against injury.



## 6. Fixed skai-padding with edge padding

Combine the fixed skai-padding with the edge padding.



## 7. Sticks- and Frame-skaipadding

Example pictures: please see page 33. All our children's care beds are available with special versions of skai-padding.



# ACCESSORIES

Equipment, which makes nursing easier




Wood care set



Easy-close-mechanism




Touch protection




Cut out




Patented lock




Patented lock II




Patented lock III




Patented lock IV




Electric adjustable side rail



270° door



Removable corner posts




Protection for the bed base

# ACCESSORIES



Children's name




Coloured sticks




Motor cover




Drawers



External lock on side of bed



Neet-roof




Feed tube hole




Feed tube hole



Entrance help



Raising aid



Play trainer wood



Play trainer metal

# Create your own children's bed



## Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

### 1. Lying area of the bed?

65 x 155 cm  70 x 140 cm  90 x 170 cm   
90 x 200 cm  100 x 200 cm  other size: \_\_\_\_\_

### 2. Which barred side do you need?

Height (cm): \_\_\_\_\_

### 3. Do you use a patient-lifter?

Yes  No

### 4. Which lying height (top of slatted base) do you need?

40 cm  50 cm  60 cm  70 cm  other height: \_\_\_\_\_

### 5. Do you need an electric height adjustment?

Yes  No

### 6. Which slatted bed base do you need?

no adjustment  manual adjustment  protection for the bed base   
electric adjustment  Trendelenburg function

### 7. What kind of doors do you want?

2-winged door  4-winged door  4-winged sliding door  270° door   
adjustable side rail  electric adjustable side rail

### 8. How many doors do you need?

front side  front and back sides  head and footelement

### 9. Locking mechanism?

patented lock I  patented lock II  patented lock III  patented lock IV

### 10. Do you need perspex?

Yes  Where? \_\_\_\_\_ No

### 11. Kind of wood?

beech  pine

### 12. Colour feature?

natural  yellow  orange  blue-purple  dark-purple   
blue  green  white stained  white lacquered  other colour: \_\_\_\_\_

coloured sticks  children's name: \_\_\_\_\_

### 13. Running gear?

wheels  without wheels

### 14. Do you want containers?

Yes  No

### 15. Do you want a mattress?

Yes  No

polyfoam mattress  interior spring mattress  cold foam mattress

Reha-Star II  Reha-Star Plus  intercontinence cover

### 16. Do you need a padding?

Yes  No

skai-padding  fixed skai-padding  skai-padding for the top edge of the cot

motif-padding: Paula  Paul  Greta

height of the padding: \_\_\_\_\_ colour of the padding: \_\_\_\_\_

### 17. Do you need Vinyl-windows?

Yes  Where? \_\_\_\_\_ No

### 18. Equipment?

touch-protection  net roof  raising aid  motor cover (Knut)

mattress-support  side up aid  removeable posts (Niklas)

- Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

- Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, data

stamp, sign



SAVI GmbH  
Reha Solutions  
Lise-Meitner-Str. 10  
D-48691 Vreden

Postfach 1123  
D-48684 Vreden  
Fon: +49 (0) 25 64 - 93 09-70  
Fax: +49 (0) 25 64 - 93 09-60

Very easy: Fax your order to  
+49 25 64 - 93 09 - 60  
or write an e-mail to  
sales@savi.de  
(Don't forget your contact details)

Or use our online-form at:  
**www.savi.de**

[www.savi.de](http://www.savi.de)  
[info@savi.de](mailto:info@savi.de)

Facebook:  
SAVI.GmbH

SAVI APP



SAVI GmbH  
Reha Solutions  
Lise-Meitner-Straße 10  
D-48691 Vreden

PO-Box 1123  
D-48684 Vreden  
Fon: +49(0)25 64 -93 09-70  
Fax: +49(0)25 64 -93 09-60

please scan QR-Code



Your dealer: