

Niklas



Be safe and free – Niklas offers your child unlimited freedom thanks to individual design options.

Our Niklas offers extreme flexibility and functionality in the barred sides. This bed can be reached from all sides because the grilles can be lowered independently of each other.

Maximum safety combined with limitless freedom – a good combination for well-being.



Electric height adjustment



Head- & footelement manual or electric adjustment (option)



More info about
Niklas:
www.savi.de/beds/niklas

EN 60601-2-52
DIN EN 32623

According to Medical Devices Act
CE-Markering
Declaration of Conformity

Niklas

Why a bed from SAVI?

- Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- Stability and reliability
- Comfortable, back-friendly height electric
- Many bed sizes
- Individual equipment options
- Optimal eye contact with your child
- Large, safe indoor area
- Colourful padding for maximum protection
- Short delivery times





Niklas' individual strengths

- All barred sides can be lowered
- 60 cm barred side height
- Removable corner posts (option)
- Low entrance (approx. 47 cm) (option)
- Trendelenburg- / Antitrendelenburg-function A+B version (option)



Choice of colours (see on page 68)



perspex + sticks ,
e.g. with adjustable sides
in deep-purple colour accent



all perspex,
e.g. with electrically adjustable sides
in deep-purple colour accent



all perspex,
e.g. with electrically adjustable sides
in deep-purple colour accent



Three kinds of our motifpadding are available.
You will find more information on our website:
www.savi.de/skai-padding

Technical data

Lying area	External dimension	Article-no.
100 x 200 cm	115 x 215 cm	131.50.006
90 x 200 cm	105 x 215 cm	131.50.006.L
90 x 170 cm	105 x 185 cm	131.50.006.M
65 x 155 cm	80 x 170 cm	131.50.006.S
Doors		
All sides adjustable in steps (lowerable) •		
Lying height (top of slatted bed base)		
76 - 116 cm electric height adjustment •		
47 - 87 cm electric height adjustment (only in combination with 99.737 possible, see page 20) o.r.		
Barred side height		Overall height
60 cm	140 cm	•
40 cm + 20 cm *	155 cm	99.737 +
Slatted bed base		
Manual adjustment	99.030	•
Electric adjustment	99.201	+
Trendelenburg- / Antitrendelenburgfunktion A (sloping position)	99.212	+
Trendelenburg- / Antitrendelenburgfunktion B (sloping position)	99.213	+
Running gear		
Four castor running gear with total locking mechanism •		
Designer wooden castors	99.885	n.a.c.
Kind of wood		
Beechwood „natural oiled“ (further colours see p. 69) •		
Mattresses		
Polyfoam mattress •		
Cold foam mattress	99.031.5	+
Reha-star mattress plus (up to stage 2)	99.031.3.PLWS	+
Reha-star mattress II (up to stage 4)	99.031.3.ZWEI	+
Incontinence cover	99.031.3ZA50	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

* Indication

The electric height adjustment 47 - 87 cm is only in combination with a barred side height of 40 cm + 20 cm high slip-possible.

Technical data

Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	n.a.c.
Orange	90.204	n.a.c.
Blue-purple	90.405	n.a.c.
Dark-purple	90.406	•
Blue	90.515	n.a.c.
Green	90.618	n.a.c.
Equipment		
Children's name on side of bed		n.a.c.
Coloured sticks		n.a.c.
PerspeX + sticks (as shown) ①	99.810LEXANNIK	+
4 removable corner posts	99.889	+
Only 1 side adjustable, 3 sides fixed	99.730	+
Only 2 sides adjustable, 2 sides fixed	99.731	+
1 electric adjustable side rail	99.735	+
2 electric adjustable side rails	99.736	+
Feed tube hole	99.840	o.r.
Raising aid	99.802	+
Sit up aid	99.820	+
Adaption of the bed for sloping ceiling		o.r.
2 drawers	99.040	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		
Motifpadding Greta, 50 cm	99.050GRETA	+
Motifpadding Paul, 50 cm	99.050PAUL	+
Motifpadding Pauline, 50 cm	99.050PAULINE	+
Skai-padding 50 cm height	99.051NIKLAS	+
Fixed skai-padding	99.051.FEST	o.r.
Fixed skai-padding with vinyl windows	99.051.FEST	o.r.

• = standard o.r. = on request + = accessories n.a.c. = no additional charge



Niklas

Key data of Niklas

- Barred side height 60 cm, all side-rails adjustable
- Trendelenburg-/Antitrendelenburg-function (option)
- Lying height 76 - 116 cm or 47 bis 87 cm (low entrance)
- Electrically adjustable sides (option)
- Removable corner posts (option)
- Perspex (optional)
- Choice of colours



perspex + sticks,
e.g. low entrance (47 - 87 cm) and two top
elements (40 + 20 cm)



Benefits of electric adjustable sides

- Easy to use at the push of a button
- Lockable control element prevents unsupervised opening
- Can also be operated by people with disabilities
- Minimal operating strength required (pushbutton)
- Infinitely adjustable



Recessed safety grip

- Easy to use manually
- The safety gap on both recesses prevents unsupervised opening
- Minimal operating strength required (slider buttons)
- Grid adjustment



recessed safety grip (basic)



Removable
corner posts

Optimum protection thanks to fixed Skai-padding

Your child is very active? With the firmly integrated Skai-padding (75 mm height), you offer your child optimum protection.

For more information on the different padding options, see page 69/70!



*Only the best
for your child!*



ID: 20841



ID: 16581



ID: 15368



ID: 17167



ID: 10352



ID: 20631

Trendelenburg-function A

Inclination of the complete bed



Trendelenburg-function B

Inclination of the slatted bed base



ID: 10394
adaption of the bed
for sloping ceiling

„I love a good
place to sleep
and lounge.“



ID: 10682
canopy bed, total height 215 cm,
entrance height 30cm

Little tip:

In this catalogue you can discover many bed variants that we have already produced. You can recognize them by the „ID“ below the illustrations.



Examples



ID: 16069
drawers



ID: 10400
coloured sticks



ID: 18692
white stained



ID: 18499
without castors



ID: 20214
closeable feed tube hole
with sliding mechanism
(reclosable)



ID: 19204
small castors



ID: 19923
drawers



Create your own SAVI®-children's bed!

Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our beds.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

- 65 x 155 cm 70 x 140 cm 90 x 170 cm
 90 x 200 cm 100 x 200 cm other size: _____

2. Which barred side do you need? Height (cm): _____

3. Do you use a patient-lifter? Yes No

4. Which lying height (top of slatted base) do you need?

- 40 cm 50cm 60 cm 70 cm other height: _____

5. Do you need an electric height adjustment? Yes No

6. Which slatted bed base do you need?

- no adjustment manual adjustment protection for the bed base
 electric adjustment Trendelenburg function

7. What kind of doors do you want?

- 4-winged door 270° door folding door
 adjustable side rail electric adjustable side rail

8. How many doors do you need?

- front side front and back sides head and footelement

9. Locking mechanism?

- SAVI-lock I SAVI-lock II SAVI-lock III SAVI-lock IV

10. Do you need perspex?

- Yes. Where?: _____ No

11. Kind of wood/wood-colour? Wood: beech pine

- Wood-colour: natural basalt wenge barrique teak

12. Colour feature?

- natural yellow orange blue-purple dark-purple blue
 green white stained white lacquered other colour: _____
 coloured sticks children's name: _____

13. Running gear? wheels without wheels

14. Do you want containers? Yes No

15. Do you want a mattress? Yes No

If yes, which one?

- polyfoam mattress interior spring mattress cold foam mattress
 Reha-Star II Reha-Star Plus intercontinence cover

16. Do you need a padding? Yes No

If yes, which one?

- skai-padding fixed skai-padding skai-padding for the top edge of the cot
 motif-padding: Pauline Paul Greta

height of the padding: _____ colour of the padding: _____

17. Do you need Vinyl-windows?

- Yes. Where? _____ No

18. Equipment?

- touch-protection net-roof raising aid motor cover (Knut)
 mattress-support side up aid feed tube hole removeable posts (Niklas)
 feed tube hole (Lexan)

Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, date

stamp, sign

Very easy:

Just send us your order via eMail
or use our online-form at: www.savi.de

SAVI GmbH Reha Solutions

Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70

D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60

