

Little Mila



Cuddly and without hurdles – with extra flat entrance of 7 cm

Due to the very low construction „Little Mila“ offers a lying height of 7 cm (up to 80 cm) including electric head and foot adjustment.



Trendelenburg- /Antitrendelenburg-
function



Barred side 70 cm



Electric height
adjustment



Head- & footelement manual or
electric adjustment (option)



More Info about
Little Mila:
www.savi.de/en/care-beds-for-children/

EN 60601-2-52
DIN EN 50637

According to Medical Device Regulation
CE-Markering
Declaration of Conformity

Little Mila

Why a bed from SAVI?

- Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- Stability and reliability
- Comfortable, back-friendly care height electric
- Many bed sizes
- Individual equipment options
- Optimal eye contact with your child
- Large, safe indoor area
- Colourful padding for maximum protection
- Short delivery times

Little Mila's individual strengths

- Barred side height 70 - 120 cm
- 8 doors
- Trendelenburg- / Antitrendelenburg-function
- Lying height of 7 cm minimum (top of slatted bed base)
- Electric head / foot adjustment
- Fixed integrated padding (option), retrofittable (just possible with Little Mila)

Lockable Handset-box (for all beds available)



Safety
first!

All model variants are available with the following barred side heights

Barred side height 70 cm,
Lying height 7 - 80 cm



Barred side height 90 cm,
Lying height 7 - 80 cm



Barred side height 120 cm,
Lying height 7 - 80 cm



Technical data

Lying area	External dimension	Article-no.
90 x 200 cm (size L)	103 x 232 cm	
Barred side height	Overall dimensions	
70 cm	83 cm	
90 cm	102 cm	
120 cm	132 cm	
Doors		
Folding doors on both sides (8 doors)		•
Locking mechanism		
SAVI-lock IV	99.034A	•
Lying height (top of slatted bed base)		
7 - 80 cm electrical		•
Slatted bed base		
Electric adjustment	99.200	•
Trendelenburg-/Antitrendelenburg-function (tilting of the entire bed))	99.213	•

• = standard o.r. = on request + = accessories

Our great fixed integrated padding
can be easily retrofitted for.

No conversion required – simply screw in – done!
(only for Little Mila)



Technical data

Running gear		
Four castor running gear made of plastic 75 mm		•
Kind of wood		
Beechwood „natural oiled“ (further colours see p. 79)		•
Mattresses		
Polyfoam mattress	99.031.4	•
Cold foam mattress	99.031.5	o.r.
Reha-star mattress plus (up to stage 2)	99.031.3.PLUS	+
Reha-star mattress II (up to stage 4)	99.031.3.ZWEI	+
Incontinence cover	99.031.3ZA50	+
Coloured features		
Natural	90.100	n.a.c.
Equipment		
Perspex		•
External lock on side of bed	99.706	+
Tube feed-through	99.840	o.r.
Raising aid	99.802	+
Sit up aid	99.820	+
Wood surface care set	99.880	+
Padding		
Motifpadding MALLU, 50 cm height	99.050MALLU	+
loose Skal-padding 50 cm height	99.051	+
Fixed Skal-padding 50 mm with integrated vinyl windows	99.051.FEST	+
Edge padding	99.053	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge



Little Mila

Optimum protection thanks to fixed Skai-padding

Your child is very active? With the fixed Skai-padding (50 mm thickness), you offer your child optimum protection.

For more information on the different padding options, see page 79/ 80!

Easy to retrofit



This Skai-padding is available in grey!



(other colours on request)



120 cm barred side height



Optimal view thanks to integrated vinyl windows

The integrated vinyl windows fit flush with the fixed Skai-padding and thus offer a similar padding effect to the permanently integrated foam padding. The windows offer an optimal view. So you keep an eye on your child easily and your child feels visibly more comfortable.



90 cm barred side height



"A bathing day in the river
is a lot of fun,
but I always get
suuuper tired."



Scan the QR code to find out more about
our fixed padding with integrated vinyl
windows.

Create your own



Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

- ☐ 65 x 155 cm ☐ 70 x 140 cm ☐ 90 x 170 cm
☐ 90 x 200 cm ☐ 100 x 200 cm other size: _____

2. Which barred side do you need? Height (cm): _____

3. Do you use a patient-lifter? ☐ Yes ☐ No

4. Which lying height (top of slatted base) do you need?

- ☐ 40 cm ☐ 50 cm ☐ 60 cm ☐ 70 cm other height: _____

5. Do you need an electric height adjustment? ☐ Yes ☐ No

6. Which slatted bed base do you need?

- ☐ no adjustment ☐ manual adjustment ☐ step-protection for the bed base
☐ electric adjustment ☐ Trendelenburg-function

7. What kind of doors do you want?

- ☐ 4-winged folding sliding door ☐ 270° door ☐ side rail "hide away" ☐ folding door
☐ lowerable side rail ☐ electric lowerable side rail

8. How many doors do you need?

- ☐ front side ☐ front and back sides ☐ head and footelement

9. Locking mechanism?

- ☐ SAVI-lock I ☐ SAVI-lock II ☐ SAVI-lock III ☐ SAVI-lock IV

10. Do you need perspex?

- ☐ Yes. Where?: _____ ☐ No

11. Kind of wood / wood-colour? Wood:

- ☐ beech ☐ pine
 Wood-colour: ☐ natural ☐ basalt ☐ wenge ☐ barrique ☐ teak ☐ walnut
☐ white stained ☐ white lacquered

12. Colour feature?

- ☐ natural ☐ yellow ☐ orange ☐ red ☐ grey ☐ blue-purple ☐ purple
☐ blue ☐ green ☐ white ☐ pink ☐ other colour: _____
☐ coloured sticks ☐ children's name: _____

13. Running gear? ☐ castors ☐ without castors

14. Do you want drawer? ☐ Yes ☐ No

15. Do you want a mattress? ☐ Yes ☐ No

If yes, which one?

- ☐ polyfoam mattress ☐ interior spring mattress ☐ cold foam mattress
☐ Reha-Star II ☐ Reha-Star Plus ☐ incontinence cover

16. Do you need a padding? ☐ Yes ☐ No

If yes, which one?

- ☐ loose Skai-padding ☐ fixed Skai-padding ☐ edge-padding ☐ motifpadding MALU
 height of the padding: _____ colour of the padding: _____

17. Do you need vinyl windows?

- ☐ Yes. Where? _____ ☐ No

18. Equipment?

- ☐ removable posts (Niklas) ☐ net-roof ☐ raising aid ☐ lockable handset box
☐ tube feed-through (gap) ☐ infusion holder ☐ drawers ☐ mattress holder
☐ touch-protection for height adjustment ☐ foot cut-out (e.g. for patient lift)
☐ central locking device for 4-roller chassis (Knut a+b only) ☐ access aid made of artificial leather
☐ supply opening (hole) with closing mechanism

☐ Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

☐ Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, date

stamp, sign

Very easy:

Just send us your order via email

or use our online-form at: www.savi.de/en/home/

SAVI GmbH Reha Solutions

Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70

D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60