

Care beds for children | Therapy couches | Travel bed | Babywalker

Welcome

Love of life for people with disabilities

I'm Günter Bußmann, and I founded this company in 1995 with the aim of manufacturing products that would provide other people with a love of life and make their everyday life easier.

The choice had to be wood as a natural product! The pleasant feel when you touch wood fascinates children in particular.

We are continuously introducing innovations into our products in collaboration with establishments such as children's homes, and in particular by talking with parents who are caring for their children.

If you have a challenge you need to solve, would like some advice on the right therapy bed for a child or want to order a bed, please do not hesitate to give us a ring:

funt. Isol

+49 (0) 25 64 / 93 09 - 70



Our aim is to reconcile quality, ergonomics and your personal wishes with nature, so that both you and your child are satisfied.



Environmentally conscious manufactory

As a manufacturer of therapy beds for children, nursing beds, heavy duty beds, Huntington's Chorea beds and therapy couches, we follow a philosophy of thinking holistically, and we also attach great value to the subject of the environment and energy when we are putting this into practice.

Materials

All the materials incorporated into our products, and especially the woods we use for producing them, come from afforestation areas and forests in Europe that are managed in an ecologically sound way. We obtain them from PEFC certified suppliers.

Production site

Our production buildings are energy-optimised and the production area is also fitted throughout with under oor heating (heat pumps) Our bright rooms are ooded with light and require hardly any additional lighting during the day, which means that we can save a lot of energy here and protect the environment!



Machinery and equipment

The refuse generated during the production process (dust and wood shavings) is compressed into briquettes and made available for burning so as to protect the environment.

The exhaust air is filtered, and by means of a heat recovery system more than 80% of the residual heat can be returned to the building. Modern, energyoptimised machinery ensures that our energy consumption is as low as possible.



Energy procurement (electricity)

A photovoltaic system of almost 200 KW is used to supply our electricity. Using this we are able to supply the energy needed for the heat pump and the majority of the energy for the production facilities.



Quality and processing

There are many factors that make a product valuable:

Raw materials

At the start of production comes the critical selection of raw materials. In our factory we only work with knotless words from an environmentally friendly perspective.

Surface

In order to maintain the natural character of the wood, the surfaces are exclusively finished with natural wax which has been tested in terms of its health and ecological aspects and in accordance with DIN EN 71 (safety of toys).

Ergonomics


When developing the beds we took care to ensure that the edges are well rounded during processing. The fittings are also on the inside to prevent injury.

Standards

All our beds comply with the European standards, bear the CE mark and are compatible with one another. We thus offer the optimal conditions for you to find the ideal solution depending on your individual requirements.

Safety

All the beds in the SAVI range offer the childproof patented lock which prevents the doors being opened inadvertently (DBGM protected).





We are specialized in nursing beds, comfort beds and special beds for seniors, heavyweights and chorea Huntington's patients.



Welcome

Advice

We find individual solutions for our customers thanks to our experience and creativity.

We offer solutions for every sized budget and will realise your desired requirements.

We would be pleased to send you a written, no obligation quotation or advise you by phone.

Tip: On page 71 you will find a customization form.

Fon: +49 (0) 25 64 / 93 09 - 70

Fax: +49 (0) 25 64 / 93 09 - 60

Custom made beds

We also produce individual custom made beds. Our technical team will be happy to talk through the possible options with you.

Showroom

Why not visit our production site in Germany? If you register in advance by phone we will schedule sufficient time to provide the optimal advice for you and/or your children. We also offer product training courses at our premises for specialist retailers!

Our clearly structured manufacturing enables us to realise many customers' wishes quickly and at a reasonable price.

Why a bed from SAVI?

- Great safety, thanks to sophisticated and safe locking systems that prevent your child from unintentionally opening from the inside
- Biologically neutral, pollutant-free wood treatment according to DIN 71-3 (safety of toys)
- The best for your child professional construction, high stability and reliability
- Comfortable, back-friendly height due to electric height adjustment (stepless)
- Many bed sizes available (up to 180 x 220 cm)
- Individual equipment options, e.g. rigid, manual or infinitely variable electric slat bed base
- Optimal eye contact with your child through the safety windows made of special perspex
- A large, safe indoor area of the bed (100 x 200 cm) offers your child enough playing space even at night and during the day and gives you the safety for your child that you need to do everyday things in peace and quiet
- Colourful padding for the interior of the bed always offers the maximum protection against injuries, even for the most active children
- Short delivery times of 2 3 weeks







content

Children's nursing beds and therapy couches from our rehabilitation production







Knut







Lola



Jonas





Felix

Travel bed | Babywalker

Contents

| SAVI Reha Solutions | p. 02 |
|---------------------|-------|
| Tom | р. 06 |
| Tom Paddy | p. 14 |
| Niklas | p. 16 |
| Knut | p. 24 |
| Knut (b) 40 | p. 36 |
| Savoir-Vivre | p. 38 |
| Jonas | p. 46 |
| Lola | p. 52 |
| Felix | p. 58 |
| Travel bed | p. 62 |
| Babywalker | p. 63 |
| Accessories | p. 64 |
| Slatted bed base | p. 66 |
| Mattresses | p. 67 |
| Range of colours | p. 68 |
| Padding | p. 69 |
| Customization form | p. 71 |

Tom



Special beds for unique children – our Tom can be designed variably according to the needs of your child and therefore always offers the right solution.

Thanks to the four-winged folding sliding door, the barred sides and the flexible height adjustment possibilities, our Tom provides an optimal therapeutic support. Due to the low entry height of a little about 20 cm even smaller children can get into the bed on their own.



low entry height of about 20 cm possible,







EN 60601-2-52 DIN EN 32623 According to Medical Devices Act CE-Markering Declaration of Conformity

Tom

Best security for your child, thanks to optimal locking systems

Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- Comfortable, back-friendly height electric
- · Many bed sizes
- · Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- · Colourful padding for maximum protection
- Short delivery times

Tom' individual strengths


- . Four-winged door to open the complete front
- Castors
- · High barred sides (option)
- . Low entry height of 20 cm minimum (optional)
- Trendelenburg-/Antitrendelenburg-function (option)
- · Fixed skai-padding (also with vinyl windows) (option)



Choice of colours (see on page 68)

A lot of safety, even for particularly active children







Barred side height 135 cm, e.g. in blue-purple colour accent with wooden castors



Technical data

| | Charles of Charles and Charles | | | 1000 |
|-------------------------|--|---|-------------------------------|--------|
| Lying area | External dimension | | Article-no. | |
| 100 x 200 cm | 115,5 x 223 cm | | 155,34,008 | |
| 90 x 200 cm | 105,5 x 223 cm | | 155.34.008.L | |
| 90 x 170 cm | 105,5 x 193 cm | | 155.34.008.M | |
| 65 x 155 cm | 80,5 x 178 cm | | 155.34.008.5 | |
| 70 x 140 cm | 85,5 x 163 cm | | 155.34.008.XS | |
| Doors | | | | |
| Four-winged do | or to open the complete | front | | |
| Four-winged do | or on both sides | | 99.750 | + |
| 270° door, door | outwardly opening | | 99.704 | + |
| 270° door, door | outwardly opening on b | ooth sides | 99.704BEID | + |
| Locking mechan | nism | | | |
| Child-proof SAV | /Flock I | | | • |
| Child-proof SAV | /Flock I metal | | 40100003 | + |
| SAVI-lock II | | | 99.033 | + |
| SAVI-lock III | | | 99.034 | + |
| SAVI-lock IV | | | 99.034A | + |
| Lying height (to | p of slatted bed base) | | | |
| | ctric height adjustment | | | |
| Barred side hei | And the second s | Overall heigh | T. | 337.0 |
| 72 cm | JANA . | 114 cm | | + |
| 113 cm | | 155 cm | | |
| 135 cm | | 177 cm | | + |
| Special sizes on | remiest | | | + |
| Slatted bed bas | | | | |
| No adjustment | *** | | | + |
| Manual adjustn | Karaka Giri | | | |
| Electric adjustm | | | 99.200 | + |
| | / Antitrendelenburg-func | tion (sloping position) | 99.211 | 4 |
| Running gear | | mana inclusio I decimaly a | | |
| | ning gear with total locki | ing mechanism | | |
| Designer woods | | | 99.885 | n,a.c. |
| EDDING MAKEN BUTCHELS & | and low entrance heigh | t of 20 cm | 99.883 | + |
| | low entrance height of | | 99.881 | + |
| Kind of wood | armorna regint or | | 22,5001 | |
| CONTRACTOR STATE | stural oiled" (further colo | ours see p. 69) | | |
| Mattresses | paramer conc | W 174 | | |
| Polyfoam mattr | ess | | | |
| Cold foam matt | | | 99.031.5 | + |
| 100 | ess plus (up to stage 2) | | 99.031.3.PLUS | + |
| | ess II (up to stage 4) | | 99.031.3.ZWEI | + |
| Incontinence co | | | 99.031.3ZA50 | + |
| | | 000000000000000000000000000000000000000 | | |
| = standard | o.r. = on request | + = accessories | n.a.c. = no additional charge | |

Technical data

| Coloured features | | |
|--|-------------|--------|
| Natural | 90.100 | n.a.c. |
| Yellow | 90.128 | n.a.c. |
| Orange | 90,204 | n.a.c. |
| Blue-purple | 90,405 | • |
| Dark-purple | 90.406 | n.a.c. |
| Blue | 90.515 | n.a.c. |
| Green | 90.618 | n.a.c. |
| Equipment | | |
| Children's name on side of bed | | n.a.c. |
| Coloured sticks | | n.a.c. |
| Perspexs instead of sticks | | • |
| Sticks instead of perspexs | 99.815 | n.a.c. |
| Touch protection for height adjustment | 99.900 | + |
| External lock on side of bed | 99.706 | + |
| Feed tube hole | 99.840 | o.r. |
| Raising aid | 99.802 | + |
| Net-roof | 99.800 | + |
| Adaption of the bed for sloping ceiling | | 0.1. |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Padding | | |
| Skai-padding 50 cm height | 99.051 | + |
| Skai-padding 100 on height | 99.051.100 | + |
| Fixed skai-padding | 99.051.FEST | O.F. |
| Fixed skai-padding with vinyl windows | 99.051.FEST | 0.1. |
| Egde padding | 99.052 | .+ |
| | | |

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Our skai-padding is in many colours available! You can find a pre-selection on our website at: www.savi.de/skai-padding

Tom



Thanks to its folding-sliding mechanism, the door offers the maximum stability and precise operation in every position. The design is a harmonious combination of quality and function. The entire front side can be opened in this way.

Tom adapts perfectly to the situation and makes everyday life with your child easier.

Naturally by request we can also offer customised additional doors, and custom-made sloping roofs are also not a problem for us.

Safety first – SAVI door locks

- Child-proof SAVI-lock I (base)
- Child-proof SAVI-lock I made of metal
- 3 SAVI-lock II
- SAVI-lock III: The labyrinth guide prevents the child opening the bed by accident. This lock is keyless and has no magnets. Can also be retrofitted.
- SAVI-lock IV: Here too the bed is opened by means of a labyrinth guide, but in addition this lock also operates against resistance. This lock is keyless, without magnets and can be retrofitted.
- External lock on side of bed (increases the stability and provides more secure)

Great solutions for door locks - highest safety first



Little tip:

In this catalogue you can discover many bed variants that we have them by the "ID" below the


Key data of Tom

- · Barred side height 72 cm, Barred side height 113 cm, Barred side height 135 cm / special sizes on request
- . Slatted bed base 34 104 cm or Slatted bed base 20 - 90 cm (without wheels)
- Trendelenburg-/Antitrendelenburg-function (option)
- · Choise between perspexs or sticks

270° door to fold back completely

- · Easy to use manually
- . Doors on the head and foot elements fold back completely
- 270° version can also be supplied with firm skai-padding





Tom

All functions of the bed can always be used without removing the padding beforehand!

Optimum protection thanks to fixed Skai-padding

Your child is very active? With the firmly integrated skai-padding (75 mm height), you offer your child optimum protection.

For more information on the different padding options, see page 69/70!

Optimal view thanks to integrated vinyl windows

The integrated vinyl windows fit flush with the firm skai-padding and thus offer a similar padding effect to the permanently integrated foam padding. The windows offer an optimal view. So you keep an eye on your child easily and your child feels visibly more comfortable.







Examples










ID: 20048









Tom Paddy

Tom Paddy's individual strengths

- . Low entry height of 20 cm minimum
- High barred sides 113 cm
- · Fixed skai-padding with vinyl windows
- · Castors (option)
- Trendelenburg-/Antitrendelenburg-function (option)







Now entry height of about 20 cm possible

Tom Paddy


Technical data

| Lying area | External dimension | Article-no. | |
|---|--|-------------------|---|
| 90 x 200 cm | 120 x 238 cm | 155.34.008.L PADD | Y |
| 90 x 170 cm | 120 x 208 cm | 155.34.008.M PADD | Y |
| Doors | | | |
| Four-winged do | or to open the complete front | | |
| Four-winged do | or on both sides | 99.750 | + |
| and the second of the second of the second of | r outwardly opening | 99.704 | + |
| 270° door, doo | r outwardly opening on both sides | 99.704BEID | + |
| Locking mechan | nism | | |
| Child-proof SA\ | | | • |
| Child-proof SAN | /Flock I metal | 40100003 | + |
| SAVI-lock II | | 99.033 | + |
| SAVI-lock III | | 99.034 | + |
| SAVI-lock IV | | 99.034A | + |
| Lying height (to | p of slatted bed base) | | |
| 20 - 90 am elec | tric height adjustment | | |
| Door height | Overall height | | |
| 123 cm | 149 cm | | + |
| Slatted bed bas | e | | |
| Manual adjustn | nent | | |
| Electric adjustm | | 99.200 | + |
| Trendelenburg- | /Antitrendelenburg-function (sloping position) | 99.211 | + |
| Running gear | | | |
| | and low entrance height of 20 cm | 99.883 | + |
| | n low entrance height of 27 cm | 99.881 | + |
| Kind of wood | | | |
| | atural oiled" (further colours see p. 69) | | • |
| Mattresses | Paris 1 20 10 | | |
| Polyfoam mattr | ess : | | • |
| Cold foam mat | | 99.031.5 | + |
| | ess plus (up to stage 2) | 99.031.3.PLUS | + |
| | ress II (up to stage 4) | 99.031.3.ZWEI | + |
| Incontinence co | wer | 99.031.3ZA50 | + |
| = standard | o.r. = on request += accessories | | |

Technical data

| Equipment | | |
|--|---------------|---|
| External lock on side of bed | 99.706 | + |
| Raising aid | 99.802 | + |
| Net-roof | 99.800 | + |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Bluetooth-control | | |
| Padding | | |
| Fixed skai-padding with vinyl windows | 99.051.FEST | |
| • = standard o.r. = on request + : | = accessories | |

Fixed padding in creme or sky-blue.



Níklas



Be safe and free - Niklas offers your child unlimited freedom thanks to individual design options.

Our Niklas offers extreme flexibility and functionality in the barred sides. This bed can be reached from all sides because the grilles can be lowered independently of each other.

Maximum safety combined with limitless freedom - a good combination for well-being.






EN 60601-2-52 DIN EN 32623 According to Medical Devices Act CE-Markering Declaration of Conformity

Níklas

Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- · Comfortable, back-friendly height electric
- Many bed sizes
- · Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- · Colourful padding for maximum protection
- · Short delivery times



Niklas' individual strengths

- · All barred sides can be lowered
- . 60 cm barred side height
- · Removable corner posts (option)
- · Low entrance (approx. 47 cm) (option)
- Trendelenburg-/Antitrendelenburg-function A+B version (option)



Choice of colours (see on page 68)




perspexs + sticks ①, e.g. with adjustable sides in deep-purple colour accent



all perspexs, e.g. with electrically adjustable sides in deep-purple colour accent



an perspexs, e.g. with electrically adjustable sides in deep-purple colour accent



Three kinds of our motifpadding are available. You will find more information on our website: www.savi.de/skai-padding

* Indication

The electric height adjustment 47 - 87 cm is only in combination with a barred side height of 40 cm + 20 cm high slip-possible.

Technical data

| Lying area | External dimension | Article-no. | | |
|-------------------|----------------------------|--|--------------|--|
| 100 x 200 cm | 115 x 215 cm | 131.50.006 | | |
| 90 x 200 cm | 105 x 215 cm | 131.50.006.L | | |
| 90 x 170 cm | 105 x 185 cm | 131.50.006.M | | |
| 65 x 155 cm | 80 x 170 cm | 131.50.006.5 | 131.50.006.5 | |
| Doors | | ACCORDED NO. | | |
| All sides adjusta | able in steps (lowerable) |) | • | |
| Lying height (to | p of slatted bed base) | | | |
| 76 - 116 cm ele | ctric height adjustment | | • | |
| 47 - 87 cm elec | ctric height adjustment(| only in combination with 99.737 possible, see page 20) | O.f. | |
| Barred side heig | ht | Overall height | | |
| 60 cm | | 140 cm | ٠ | |
| 40 cm + 20 cm | • | 155 cm 99.737 | + | |
| Slatted bed base | | | | |
| Manual adjustm | ent | 99.030 | ٠ | |
| Electric adjustm | ent | 99.201 | + | |
| Trendelenburg- | Antitrendelenburgfunk | tion A (sloping position) 99.212 | + | |
| Trendelenburg- | Antitrendelenburgfunk | tion B (sloping position) 99.213 | + | |
| Running gear | | | | |
| Four castor runn | ning gear with total lock | ring mechanism | | |
| Designer woode | en castors | 99.885 | n.a.c | |
| Kind of wood | | | | |
| Beechwood "na | itural oiled" (further col | ours see p. 69) | | |
| Mattresses | | | | |
| Polyfoam mattre | 255 | | • | |
| Cold foam matt | 97.5 | 99.031.5 | + | |
| Reha-star mattr | ess plus (up to stage 2) | 99.031.3.PLUS | + | |
| Reha-star mattr | ess II (up to stage 4) | 99.031.3 ZWE | + | |
| Incontinence co | ver | 99.031.3ZA50 | + | |
| = standard | | + = accessories n.a.c. = no additional charge | | |

Technical data

| Coloured features | | |
|--|----------------|--------|
| Natural | 90.100 | n.a.c. |
| Yellow | 90.128 | n.a.c. |
| Orange | 90,204 | n.a.c. |
| Blue-purple | 90,405 | n.a.c. |
| Dark-purple | 90.406 | • |
| Blue | 90.515 | n.a.c. |
| Green | 90.618 | n.a.c. |
| Equipment | | |
| Children's name on side of bed | | n.a.c. |
| Coloured sticks | | n.a.c. |
| Perspexs + sticks (as shown) ① | 99.810LEXANNIK | + |
| 4 removable comer posts | 99.889 | + |
| Only 1 side adjustable, 3 sides fixed | 99.730 | + |
| Only 2 sides adjustable, 2 sides fixed | 99.731 | + |
| 1 electric adjustable side rail | 99.735 | + |
| 2 electric adjustable side rails | 99.736 | + |
| Feed tube hole | 99.840 | 0.1. |
| Raising aid | 99.802 | + |
| Sit up aid | 99.820 | + |
| Adaption of the bed for sloping ceiling | | о.г. |
| 2 drawers | 99.040 | + |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Padding. | | |
| Motifpadding Greta, 50 cm | 99.050GRETA | + |
| Motifpadding Paul, 50 cm | 99.050PAUL | + |
| Motifpadding Pauline, 50 cm | 99.050PAULINE | + |
| Skai-padding 50 cm height | 99.051 NIKLAS | + |
| Fixed skai-padding | 99.051.FEST | o.r. |
| Fixed skai-padding with vinyl windows | 99.051.FEST | 0.1. |
| | | |





Níklas

Key data of Niklas

- · Barred side height 60 cm, all side-rails adjustable
- Trendelenburg-/Antitrendelenburg-function (option)
- Lying height 76 116 cm or 47 bis 87 cm (low entrance)
- · Electrically adjustable sides (option)
- · Removable corner posts (option)
- · Perspexs (optional)
- · Choice of colours



Infinitely adjustable



Removable

corner posts

Recessed safety grip

- · Easy to use manually
- The safety gap on both recesses prevents unsupervised opening
- Minimal operating strength required (slider buttons)
- Grid adjustment



recessed safety grip (basic)



Benefits of electric adjustable sides

· Lockable control element prevents unsupervised opening

· Can also be operated by people with disabilities

. Easy to use at the push of a button

Optimum protection thanks to fixed Skai-padding

Your child is very active? With the firmly integrated Skai-padding (75 mm height), you offer your child optimum protection.

For more information on the different padding options, see page 69/70!





Only the best for your child!











Níklas

Trendelenburg-function A

Inclination of the complete bed



Trendelenburg-function B

Inclination of the slatted bed base





Examples







Knut



cuddle safely in Knut – the bed with a high degree of safety.

Thanks to the four-winged folding sliding door, the barred sides and the flexible height and inclination adjustment possibilities, this bed provides very good therapeutic support.



Trendelenburg-/Antitrendelenburgfunction









EN 60601-2-52 DIN EN 32623 According to Medical Devices Act CE-Markering Declaration of Conformity

Knut

Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- · Comfortable, back-friendly height electric
- · Many bed sizes
- · Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- · Colourful padding for maximum protection
- · Short delivery times



- Barred side 64 89 cm
- Trendelenburg-/Antitrendelenburg-function
- Drawer/stowage (option)
- . Doors or adjustable side rails (option)
- . Low entry height of 21 cm minimum (option)



All models are available with alternative side rail heights

Knut (a)

Barred side height 64 cm lying height 41 - 76 cm or 31 - 66 cm (small castors)

Knut (b)

Barred side height 78 cm lying height 41 - 91 cm or 31 - 81 cm (small castors)

Knut (c)

Barred side height 89 cm lying height 31 - 101 cm or 21 - 91 cm (small castors)

* Indication

With adjustable side rails the lying height changes to 41 - 111cm!







Technical data

| Lying area | | External dimens | sion | Article-no. | |
|---------------------------|--|--|---|---------------|-----------|
| 100 x 200 cm | | 113 x 218 cm | | | |
| 90 x 200 cm (size L) | | 103 x 218 cm | | | |
| 90 x 170 cm (size M) | | 103 x 188 cm | | | |
| 70 x 140 cm (size XS) | | 83 x 158 cm | (Knut (c) | 103 x 158 cm) | |
| 65 x 155 cm (size 5) | | 78 x 173 cm | (Knut (c) | 98 x 173 cm) | |
| Doors | | -24-000-00-00-00-00-00-00-00-00-00-00-00-0 | AND SECURITION OF THE PARTY OF | | |
| 4-winged door to open th | | | | | • |
| Four-winged door on both | n sides | | | 99.750 | + |
| 270° door, door outwardi | Appeal of the contract of the contract of the behalf of the contract of the co | | | 99.704 | + |
| 270° door, door outward | y opening on both sides | | | 99.704BEID | + |
| Locking mechanism | | | | | |
| Child-proof SAVI-lock I | | | | | • |
| Child-proof SAVI-lock I m | etal | 40100003 | | | + |
| SAVI-lock II | | 99.033 | | | + |
| SAVI-lock III | | 99.034 | | | + |
| SAVI-lock IV | | 99.034A | | | + |
| Lying height | | | | | |
| Top of slatted bed base | | Barred side height | | | |
| (a) 41- 76 cm electric | | 64 cm (2) | | | • Knut (a |
| (b) 41- 91 cm electric | | 78 cm ⁽²⁾ | | | • Knut (b |
| (c) 31-101 cm electric | | 89 cm (1) | | | • Knut (c |
| Additional charge for 1 a | djustable side-rail | Bar height | | | |
| (+14 cm external lenght) | manual | 64 cm | | 99.710 | + |
| | electric | | | 99.710EL | + |
| | manual | 78 cm | | 99.712 | + |
| | electric | | | 99.712EL | + |
| | manual | 89 cm | | 99.714 | + |
| | electric | | | 99.714EL | + |
| Additional charge for 2 a | | Bar height | | | |
| (+14 cm external lenght) | manual | 64 cm | | 99.711 | + |
| | electric | | | 99.711EL | + |
| | manual | 78 cm | | 99.713 | + |
| | electric | | | 99.713EL | + |
| | manual | 89 cm | | 99.715 | + |
| | electric | | | 99.715EL | + |
| Slatted bed base | | | | | |
| Manual adjustment | | | | 99.030 | • |
| Electric adjustment | | | | 99.200 | + |
| Trendelenburg-/Antitrend | lelenburg-function (sloping p | osition) | | | • |
| = standard o.r. = e | on request + = accesso | ories | | | |

(1) Indication

In combination with adjustable side rails not possible.

(2) Indication

All models are available with alternative side rail heights. We advise you gladly.

Technical data

| Running gear | | |
|---|---------------|-------|
| Four castor running gear with total locking mechanism | | |
| Designer wooden castors | 99.885 | n.a.c |
| Without or small rolls ⁽¹⁾ , (21-91 cm) | 99.882 | n.a.c |
| Kind of wood | | |
| Beechwood "natural oiled" (further colours see p. 69) | | • |
| Mattresses: | | |
| Polyfoam mattress | | |
| Cold foam mattress | 99.031.5 | 0.1. |
| Reha-star mattress plus (up to stage 2) | 99.031.3.PLUS | + |
| Reha-star mattress II (up to stage 4) | 99.031.3.ZWEI | + |
| Incontinence cover | 99.031.3ZA50 | + |
| Coloured features | | |
| Natural | 90.100 | n.a.c |
| Yellow | 90.128 | n.a.c |
| Orange | 90.204 | • |
| Blue-purple | 90.405 | n.a.c |
| Dark-purple | 90.406 | n.a.c |
| Blue | 90.515 | n.a.c |
| Green | 90.618 | n.a.c |
| Equipment | | |
| Children's name on side of bed | | n.a.c |
| Coloured sticks | | n.a.c |
| Perspexs instead of sticks | | 0.1. |
| Perspexs in the lower area of the head / foot end (standard on C) | | 0.1. |
| External lock on side of bed | 99.706 | + |
| Feed tube hole | 99.840 | 0.1. |
| Raising aid | 99.802 | + |
| Sit up aid | 99.820 | + |
| 2 drawers | 99.040 | + |
| Motor cover | 99.803 | + |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Padding | a resident | |
| Motifpadding Greta, 50 cm | 99.050GRETA | + |
| Motifpadding Paul, 50 cm | 99.050PAUL | + |
| Motifpadding Pauline, 50 cm | 99.050PAULINE | + |
| Skai-padding 50 cm height | 99.051 | + |
| Fixed skai-padding | 99.051.FEST | 0.1. |
| Fixed skai-padding with vinyl windows | 99.051.FEST | 0.1. |
| Egde padding | 99.052 | + |

*Indication

If you choose KNUT C with drop down sides (artikel-no.: 99.710 - 99.715) the lying height will change to 41 - 111 cm.

Knut

Open side front thanks to four-winged sliding door

Thanks to its folding-sliding mechanism, the door offers the maximum stability and precise operation in every position. The entire front side can be opened in this way. The design is a harmonious combination of quality and function. Naturally by request we can also offer customised additional doors, and custom-made sloping roofs are also not a problem for us.







Key data of Knut

- Barred side height 64 cm,
 Barred side height 78 cm,
 Barred side height 89 cm,
 in all model variants or different height, on request
- · Lying height

Knut (a)

41 - 76 cm (with small castors 31 - 66 cm)

Knut (b)

41 bis 91 cm (with small castors 31 - 81 cm)

Knut (c)

31 bis 101 cm (with small castors 21 - 91 cm)


Hinweis: With adjustable side rails the lying height changes to 41 cm!

- Trendelenburg-/Antitrendelenburg-function
- Perspexs (option)
- · Four-winged sliding door
- Adjustable side rail (option)










Tubing furnishing

Our tubing furnishing is available as follows:







tubing on the middle foot element





ID: 17139





Knut

Adjustable side rails







Recessed safety grip

- Easy to use manually
- The safety gap on both recesses prevents unsupervised opening
- Minimal operating strength required (slider buttons)
- Grid adjustment









Versions adjustable side rails

Our adjustable side rail is available as follows:



Version 1

Adjustable side rail with 2 top elements (choise manually or electrically adjustable).



Version 2

Adjustable side rail with 1 top element. You can remove the top element if you drop down the side rail.



Special version 3

Adjustable side rail with hinged top element. Attention: This version is only possible in combination with an electric side rail!



- Easy to use at the push of a button
- · Lockable control element prevents unsupervised opening
- · Can also be operated by people with disabilities
- · Minimal operating strength required (pushbutton)
- Infinitely adjustable





Knut

270° door to fold back completely

- Easy to use manually
- . Doors on the head and foot elements fold back completely
- 270° version can also be supplied with firm skai-padding



Examples

Indication

Our children's therapy beds come in a choice of white painted or white stained (wood grain visible) finishes.









Knut

Optimum protection thanks to fixed skai-padding

Your child is very active? With the firmly integrated skai-padding (75 mm height), you offer your child optimum protection.

For more information on the different padding options, see page 69/70!



Our skai-padding is in many colours available! You can find a pre-selection on our website at: www.savi.de/skai-padding











Optimal view thanks to integrated vinyl windows

The integrated vinyl windows fit flush with the firm skai-padding and thus offer a similar padding effect to the permanently integrated foam padding. The windows offer an optimal view. So you keep an eye on your child easily and your child feels visibly more comfortable.



Tubing furnishing

On request, we manifacture all care beds, with tubing furnishing.




Edge padding



Stick-/Framepadding







Knut (b) 40

Technical data

| Lying area | External dimension | |
|--|--------------------|-----|
| 90 x 200 cm (size L) | 106,5 x 232 cm | |
| 90 x 170 cm (size M) | 106,5 x 202 cm | |
| Lying height Top of slatted bed base | Barred side height | |
| 41 - 91 am electric | 64 cm bzw. 40 cm | |
| Adjustable side rail | Bar height | |
| (+14 cm external lenght) manual | 64 cm 99.710 | |
| electric | 64 cm 99.710EL | - + |
| 2 adjustable side rails | Bar height | |
| (+14 cm external lenght) manual | 64 cm 99.711 | + |
| electric | 64 cm 99.711EL | + |
| Slatted bed base | | |
| Manual adjustment | 99.030 | |
| Electric adjustment | 99.200 | + |
| Trandalarhura / Antitrandalarhura function | /cloning position) | |

• = standard o.r. = on request + = accessories



Technical data

| Beechwood "natural oiled" (further colours see p. 69) | | • |
|---|---|--------|
| Mattresses | | |
| Polyfoam mattress | | • |
| Cold foam mattress | 99.031.5 | + |
| Reha-star mattress plus (up to stage 2) | 99.031.3.PLUS | + |
| Reha-star mattress II (up to stage 4) | 99.031.3.ZWEI | + |
| Incontinence cover | 99.031.3ZA50 | + |
| Coloured features | | |
| Natural | 90.100 | n.a.c. |
| Yellow | 90.128 | n.a.c. |
| Orange | 90.204 | • |
| Blue-purple | 90.405 | n.a.c. |
| Dark-purple | 90.406 | n.a.c. |
| Blue | 90.515 | n.a.c. |
| Green | 90.618 | n.a.c. |
| Equipment | | |
| Children's name on side of bed | | n.a.c. |
| Coloured sticks | | n.a.c. |
| Perspexs instead of sticks | | 0.1. |
| Perspexs in the lower area of the head / foot end (standard on C) | | 0.1. |
| Feed tube hole | 99.840 | 0.1. |
| Raising aid | 99.802 | + |
| Sit up aid | 99.820 | + |
| 2 drawers | 99.040 | + |
| Motor cover | 99.803 | + |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Bluetooth control | *************************************** | |
| Padding Padding | | |
| Motifpadding Greta, 50 cm | 99.050GRETA | + |
| Motifpadding Paul, 50 cm | 99.050PAUL | + |
| Motifpadding Pauline, 50 cm | 99.050PAULINE | + |
| Skai-padding 50 cm height | 99.051 | + |
| Fixed skai-padding | 99.051.FEST | 0.1. |
| Fixed skai-padding with vinyl windows | 99.051.FEST | O. F. |
| | 99.052 | |

Savoir-Vivre



Strong as an animal through everyday life - our ergonomic beds offer your child the ideal conditions for a healthy sleep.

This model combines safety with comfort. The electric lying-height adjustment is especially useful when it comes to the daily care of handicapped people.









EN 60601-2-52 DIN EN 32623 According to Medical Devices Act CE-Markering Declaration of Conformity

Savoir-Vivre


Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- · Comfortable, back-friendly height electric
- · Many bed sizes
- · Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- · Colourful padding for maximum protection
- · Short delivery times


Savoir-Vivre's individual strengths

- · Barred sides 60 220 cm (optional)
- Trendelenburg-/Antitrendelenburg-function version B (option)
- Drawer/stowage (option)











e.g. in blue-purple colour accent









Technical data

| | | OTOMORIO NO PORTO DO COMO DE C | THE PERSON NAMED IN COLUMN NAM | |
|--|--|--|--|--------|
| Lying area | | xternal dimension | | |
| 100 x 200 cm | | 13 x 206 cm | please specify size | |
| 90 x 200 cm (si | 10 Page 10 Pag | 03 x 206 cm | please specify size | |
| 90 x 170 cm (si | | 03 x 176 cm | please specify size | |
| 65 x 155 cm (si | re 5) | 78 x 161 cm | please specify size | |
| Doors | | | | |
| | open the complete front | | | • |
| Four-winged door | | | 99.750 | + |
| | outwardly opening | | 99.704 | + |
| A CONTRACTOR OF THE PARTY OF TH | outwardly opening on both sides | | 99.704BEID | + |
| Locking mechanis | | | | |
| Child-proof SAVI | | | | • |
| Child-proof SAVI- | lock I metal | | 40100003 | + |
| SAVI-lock II | | | 99.033 | + |
| SAVI-lock III | | | 99.034 | + |
| SAVI-lock IV | | | 99.034A | + |
| Lying height (top | of slatted bed base) | | | |
| 40 - 80 cm electr | ic height adjustment | | | • |
| Barred side heigh | t Overall height | | | |
| 60/70 cm | 110 cm | | 110.50.504.99101 | |
| 100/110 cm | 150 cm | | 150.50.504.99101 | |
| | for size 90 x 170 cm | | | |
| 125/135 cm | 175 cm | | 175.50.504.99101 | |
| | for size 90 x 170 cm | | | |
| 150/160 cm | 200 cm | | 200.50.504.99101 | |
| 170/180 cm | 220 cm | | 220.50.504.99101 | |
| Special sizes on re | equest | | | |
| Slatted bed base | | | | |
| Manual adjustme | nt | | 99.030 | • |
| Electric adjustme | nt | | 99.200 | + |
| Trendelenburg-7/ | Intitrendelenburgfunktion B (Schrägstellun | ig der Liegefläche | 99.213 | + |
| Running gear | | | | |
| Four castor runni | ng gear with total locking mechanism | | 99.887 | + |
| Designer wooden | castors | | 99.886 | + |
| Cut out on the fro | ont side (notify height please) | | 99.830 | n.a.c. |
| Kind of wood | | | | |
| Pinewood | | | | n.a.c. |
| Beechwood "nat | ural oiled" (further colours see p. 69) | | | |
| Mattresses | | | | |
| Polyfoam mattres | 5 | | | |
| Cold foam mattre | SS | | 99.031.5 | + |
| Reha-star mattres | s plus (up to stage 2) | | 99.031.3.PLUS | + |
| Reha-star mattres | is II (up to stage 4) | | 99.031.3.ZWEI | + |
| Incontinence cove | | | 99.031.3ZA50 | + |
| = standard | o.r. = on request += accessories | n.a.c. = no | additional charge | |

Technical data

| Coloured features | | |
|--|---------------|--------|
| Natural | 90.100 | n.a.c. |
| Yellow | 90.128 | n.a.c. |
| Orange | 90.204 | n.a.c. |
| Blue-purple | 90.405 | n.a.c. |
| Dark-purple | 90.406 | n.a.c. |
| Blue | 90.515 | n.a.c. |
| Green | 90.618 | • |
| Equipment | | |
| Children's name on side of bed | | n.a.c. |
| Coloured sticks | | n.a.c. |
| Perspexs instead of sticks | | 0.1. |
| Feed tube hole | 99.840 | 0.1. |
| External lock on side of bed | 99.706 | + |
| Raising aid | 99.802 | + |
| Sit up aid | 99.820 | + |
| Net-roof | 99.800 | + |
| Adaption of the bed for sloping ceiling | | 0.1. |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Padding : | | |
| Motifpadding Greta, 50 cm | 99.050GRETA | + |
| Motifpadding Paul, 50 cm | 99.050PAUL | + |
| Motifpadding Pauline, 50 cm | 99.050PAULINE | + |
| Skai-padding 50 cm height | 99.051 | + |
| Skai-padding 100 cm height | 99.051.100 | + |
| Fixed skai-padding | 99.051.FEST | 0.1. |
| Fixed skai-padding with vinyl windows | 99.051.FEST | 0.1. |
| Egde padding | 99.052 | + |
| | | |

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Our skai-padding is in many colours available! You can find a pre-selection on our website at: www.savi.de/skai-padding



Savoir-Vivre

Key data of Savoir-Vivre

- Barred side height 60 170 cm, other heights on request
- Electric lying height 40 80 cm
- · Four-winged sliding door
- · Four-castors running gear (option)
- Child-proof SAVI-lock
- Perspexs (optiona)



• 135 cm

• 160 cm

• 180 cm

Total heights: 110 cm | 150 cm | 175 cm | 200 cm | 220 cm

Heights

- Barred side Door
- 60 cm 70 cm • 100 cm • 110 cm
- 125 cm
- 150 cm
- 170 cm

- Total
- 110 cm
- 150 cm
- 175 cm
- 200 cm
- 220 cm



Locking of doors



external lock on side of bed (option)

heavy, massive, good



Cut out



Height: 8 cm (other heights are possible, e.g. for use of patient lifts)

great for the back

Examples











Savoir-Vivre

Optimum protection thanks to fixed skai-padding

Your child is very active? With the firmly integrated skai-padding (75 mm height), you offer your child optimum protection.

For more information on the different padding options, see page 69/70!







"Sometimes it feels so good to

stretch my hooves properly."









Jonas



On a dream journey with Jonas – our handmade beds offer your child an unmistakable sleeping comfort and particularly beautiful dreams.

This bed offers a high degree of safety. Guaranteed by the high barred sides and our proven SAVI-lock.





EN 60601-2-52 DIN EN 32623 According to Medical Devices Act CE-Markering Declaration of Conformity

Jonas

Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- · Comfortable, back-friendly height electric
- Many bed sizes
- · Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- · Colourful padding for maximum protection
- · Short delivery times



ID: 11071 Barred side height 110 cm e.g. with drawers, in blue colour accent

Jonas' individual strengths

- · Cost-effective, as no height adjustment is required
- . Barred side 60 220 cm (option)
- Drawer/stowage (option)



Technical data

| Lying area | External dimension | Article-no. | |
|---|--------------------------|------------------------|-------|
| 100 x 200 cm | 113 x 206 cm | | |
| 90 x 200 cm | 103 x 206 cm | | |
| Doors | | | |
| 4-winged door to open the complete f | ront | | |
| Four-winged door on both sides | | 99.750 | + |
| 270° door, door outwardly opening | | 99.704 | + |
| 270° door, door outwardly opening or | both sides | 99.704BEID | + |
| tocking mechanism | | | |
| Child-proof SAVI-lock I | | | |
| Child-proof SAVI-lock I metal | | 40100003 | + |
| SAVI-lock II | | 99.033 | + |
| SAVI-lock III | | 99.034 | + |
| SAVI-lock IV | | 99.034A | + |
| tying height (top of slatted hed base) | | | |
| 50 cm | | | |
| Barred side height | Overall height | | |
| 100 cm | 150 cm | 150.50.004 | |
| 125 cm | 175 cm | 175.50.004 | |
| Special sizes on request | | | |
| Slatted bed base | | | |
| Manual adjustment | | 99.030 | + |
| Electric adjustment (keine Höhenverste | ellung) | 99.200 | + |
| Running gear | | | |
| Four castor running gear with total loc | king mechanism | | 0.1. |
| Designer wooden castors | | | 0.f. |
| Lind of wood | | | |
| Pinewood | | | n.a.c |
| Beechwood "natural oiled" (further co | olours see p. 69) | | |
| | | | |
| Polyfoam mattress | | 99.031.4 | • |
| Cold foam mattress | | 99.031.5 | + |
| Reha-star mattress plus (up to stage 2 |) | 99.031.3.PLUS | + |
| Reha-star mattress II (up to stage 4) | | 99.031.3.ZWEI | + |
| Incontinence cover | | 99.031.3ZA50 | + |
| = standard o.r. = on request | + = accessories n.a.c. : | = no additional charge | |





Technical data

| Coloured features | | |
|---|---------------|--------|
| Natural | 90.100 | n.a.c. |
| Yellow | 90.128 | n.a.c. |
| Orange | 90.204 | n.a.c. |
| Blue-purple | 90.405 | n.a.c. |
| Dark-purple | 90.406 | n.a.c. |
| Blue | 90.515 | |
| Green | 90.618 | n.a.c. |
| | | |
| Children's name on side of bed | | n.a.c. |
| Coloured sticks to go to the coloured feature | | n.a.c. |
| Perspexs instead of sticks | | 0.1. |
| External lock on side of bed | 99.706 | + |
| Feed tube hole | 99.840 | O.F. |
| Raising aid | 99.802 | + |
| Adaption of the bed for sloping ceiling | | o.r. |
| 2 drawers | 99.040 | + |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Padding | | |
| Motifpadding Greta, 50 cm | 99.050GRETA | + |
| Motifpadding Paul, 50 cm | 99.050PAUL | + |
| Motifpadding Pauline, 50 cm | 99.050PAULINE | + |
| Skai-padding 50 cm height | 99.051 | + |
| Skai-padding 100 on height | 99.051.100 | + |
| Fixed skai-padding | 99.051.FEST | о.г. |
| Fixed skai-padding with vinyl windows | 99.051.FEST | o.r. |
| Egde padding | 99.052 | + |
| | | |

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Our skai-padding is in many colours available! You can find a pre-selection on our website at: www.savi.de/skai-padding



Jonas

Key data of Jonas

- Barred side height 100 125 cm, or even highter (i.e. 220 cm) available
- · Fixed lying height 50 cm (entrance height)
- · Drawers (optional)
- Child-proof SAVI-lock
- · Perspexs (option)





ID: 18397 Lower section perspexs, and top section sticks







Examples





Small heroes very big - Lola is the ideal bed for toddlers.

This bed is ideal for small children because it has a fixed, high basic lying height, making nursing more comfortable.

Its also easy on the back, thereby seeing to restful sleep.





EN 60601-2-52 DIN EN 32623 According to Medical Devices Act CE-Markering Declaration of Conformity

Lola

Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- · Comfortable, back-friendly height electric
- Many bed sizes
- · Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- · Colourful padding for maximum protection
- · Short delivery times

Lola's individual strengths

- · Good working height/slatted bed base 85 cm
- · Cupboards (option)



Choice of colours (see on page 68)







barred side height 65 cm and cupboards e.g. in green colour accent

Technical data


| Lying area | External dimension | Article-no. | |
|--|--------------------------|---|-------|
| 100 x 200 cm | 113 x 206 cm | | |
| 90 x 200 cm | 103 x 206 cm | | |
| Doors | | | |
| 4-winged door to open the complete | front | | • |
| Four-winged door on both sides | | 99.750 | + |
| 270° door, door outwardly opening | | 99.704 | + |
| 270° door, door outwardly opening of | on both sides | 99.704BEID | + |
| Locking methanism | | | |
| Child-proof SAVI-lock I | | | • |
| Child-proof SAVI-lock I metal | | 40100003 | + |
| SAVI-lock II | | 99.033 | + |
| SAVI-lock III | | 99.034 | + |
| SAVI-lock IV | | 99.034A | + |
| Lying height (top of slatted bed base |) | | |
| 85 cm | | | • |
| Barred side height | Overall height | | |
| 65 cm | 150 cm | 150.85.002 | • |
| 90 cm | 175 cm | 175.85.002 | + |
| Special sizes on request | | | + |
| Slatted bed base | | | |
| Manual adjustment | | 99.030 | + |
| Electric adjustment | | 99.200 | + |
| Running gear | Trees and the second | | |
| Four castor running gear with total lo | ocking mechanism | | 7.0 |
| Designer wooden castors | | | 0.f. |
| Kind of wood | | | |
| Pinewood | | | n.a.c |
| Beechwood "natural oiled" (further | colours see p. 69) | | |
| Mattresses | | 20.000000000000000000000000000000000000 | |
| Polyfoam mattress | | 99.031.4 | • |
| Cold foam mattress | | 99.031.5 | + |
| Reha-Star Plus (Stadium 2) | | 99.031.3 PLUS | + |
| Reha-Star II (Stadium4) | | 99.031.3.ZWEI | + |
| Incontinence cover | | 99.031.3.ZA50 | + |
| • = standard o.r. = on request | + = accessories n.a.c. = | no additional charge | |


Technical data

| Coloured features | | |
|---|---------------|--------|
| Natural | 90.100 | n.a.c. |
| Yellow | 90.128 | • |
| Orange | 90.204 | n.a.c. |
| Blue-purple | 90.405 | n.a.c. |
| Dark-purple | 90.406 | n.a.c. |
| Blue | 90.515 | n.a.c. |
| Green | 90.618 | n.a.c. |
| Equipment | | |
| Children's name on side of bed | | n.a.c. |
| Coloured sticks to go to the coloured feature | | n.a.c. |
| Perspexs instead of sticks | | 0.1. |
| External lock on side of bed | 99.706 | + |
| Feed tube hole | 99.840 | o.r. |
| Raising aid | 99.802 | + |
| Sit up aid | 99.820 | + |
| Adaption of the bed for sloping ceiling | | 0.1. |
| Cupboards (2 pieces) | 99.041 | + |
| Mattress holder | 99.891 | + |
| Protection for the bed base in case of kicks | 99.892 | + |
| Wood care set | 99.880 | + |
| Padding | | |
| Motifpadding Greta, 50 cm | 99.050GRETA | + |
| Motifpadding Paul, 50 cm | 99.050PAUL | + |
| Motifpadding Pauline, 50 cm | 99.050PAULINE | + |
| Skai-padding 50 cm height | 99.051 | + |
| Fixed skai-padding | 99.051.FEST | 0.1. |
| Fixed skai-padding with vinyl windows | 99.051.FEST | 0.1.0 |
| Egde padding | 99.052 | + |

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Lola







As flexible as Felix - the therapeutig couch without limits

This therapeutic couch offers your child a large therapeutic area and a low entry height. This makes the therapy particularly pleasant for parents, carers and children. It also conceals sufficient storage space for therapy needs and much more.

Its homely design and disinfectability makes Felix thoroughly suitable for everyday use.



Felix

Felix' individual strengths

It is often the small practical aids that decide whether something is good or bad in everyday use. As manufacturers of furniture for rehabilitation purposes, we have incorporated this valuable experience in the design of our new therapeutic bed. The result is a multifunctional couch in a design that appeals to children.

Handy: storage box for utensils like, for example, blankets and preparations

Materials: Powder-coated tubular steel frame, padded lying surface: Skai cover

Technical data - Felix (a)

Lying area (min. 130 x 100 cm, max. 218 x 100 cm)

Without castors, height adjustment 68 - 103 cm

Castors, height adjustment 75 - 110 cm

9910.75110







Technical data - Felix (b)

| Therapeutic couch (max. load 150 kg) | | |
|--|-----------------|--------|
| Lying area (min. 172,5 x 100 cm, max. 261 x 100 cm) | Article-no. | |
| Without castors, height adjustment 45 - 86 cm | 9911.4586 | |
| Castors, height adjustment 52 - 93 cm | 9911.5293 | |
| Therapeutic couch no folding sides (only fixed lying area) | | |
| Lying area | Article-ng. | |
| 70 x 140 cm | 9911.529370140 | |
| 90 x 170 cm | 9911.529390170 | |
| 90 x 200 cm | 9911.529390200 | |
| 100 x 200 cm | 9911.5293100200 | |
| Without castors, height adjustment of 45 - 86 cm | | n.a.c. |
| Castors, height adjustment of 57 - 93,5 cm | | |
| | | |

• = standard n.a.c. = no additional charge



Travel bed

Why you should purchase our travel bed

Our travel bed is made of beech wood. Quick assembly with fewer screws simplifies assembly. Our travel bed consists of:

- 1 beechwood frame
- 1 net-element, very tear-resistant
- · anti-climb protection integrated
- 1 bag

Remarkable features of our travel bed

- Dismantled size: 120 x 15 x 20 cm (l x w x h), packed in a bag without mattress
- Mounted size: 196 x 76 x 120 cm (l x w x h)
- Total weight: 13 kg

Technical data

| Lying area | Barred side height | Article-no. | |
|------------------------------|--------------------|---------------|---|
| SAVI-travel bed 70 x 190 cm | 120 am | 120,070,01,05 | |
| Mattress 70 x 190 cm, 3-fach | | | + |
| Net in black | | | |
| Net in blue | | | |
| Net in pink | | | |

• = standard + = accessories

Practically packed

in a bag





Babywalker

Why you should purchase our babywalker

Super for learning to walk! Load into, unload, transport! The child is able to transport himself!



Remarkable features of our Babywalker







Height adjustmentable hand grip

Light castors

Sharpening brake

Technical data

| Height | Width | Length | | |
|------------|---------------|-----------------------|-----|--|
| 84 cm | 54 cm | 70 cm | 3.0 | |
| Running o | gear | | | |
| 4 Light ca | stors with sh | arpening brake | (i) | |
| Storage b | OX | | | |
| Large load | ding area | | • | |
| Hand grip | | TAX AV ROOM | | |
| Height ad | justment, so | that the walker grows | 3● | |
| Load | 1 10 | 950 | | |
| Up to 30 | kg | _ | | |
| Surface | | | | |
| According | to DIN 71-3 | (security for toys) | • | |

• = standard + = accessories



Accessories

Equipment, which makes nursing easier



Wood care set



SAVI-lock I



Electric adjustable side rail




Easy-close-mechanism



SAVI-lock II



270° door



Touch protection





Removable corner posts



Cut out



SAVI-lock IV



Protection for the bed base



Accessories



Children's name



External lock on side of bed



Entrance help



Coloured sticks



Neet-roof



Raising aid



Motor cover



Feed tube hole with sliding mechanism



Play trainer wood



Drawers



Feed tube hole



Slatted bed base

Slatted bed base no adjustment (fixed):

Slats frame with individual central zone reinforcement.



Electric slatted bed base (no height adjustment):

Continuously variable hardness regulation with electric adjustment of the head, torso, pelvic, knee and foot



Fixed slatted bed base:

Base plate made of plywood (no bed frames).



Penetration protection wood

Installation above the bar frame to the slats.



very simple and clear



Mounted below the slatted frame. Spring action and reclining comfort are maintained.



Mattresses

Polyfoam mattress

Material: Pu foam

Strength: Density (35 kg/m³)

Reference substance: Blue elastic cover in cotton fabric (at

60 °C washable), optionally also in

difficult flammable execution.

Overall height

of the mattress: ca. 12 cm

Cold foam mattress

Height: Approx. 14 cm

Coating: Double jersey (65% polyester, 35% viscose)
Padding: On both sides ca. 350 g/m² climate fiber

(100% polyester), fleece

Care: Cover removable by 4-sided zipper. Halves divisible

and individually washable up to 60°C in gentle cycle with mild detergent, suitable for tumble

dryers.

Core: 7-zones-Ortho-cel®-cold-foam-core

(height of core approx. 12 cm)

Incontinence cover

Fabric: 100% polyester Coating: 100% polyurethan


Supple, skin-friendly, anti-allergic, bacteria-proof,


resistant to blood, urine, oil, grease, disinfectants

Care: Machine wash 95 °C suitable









Reha-Star mattress plus up to stage 2 (nach Prof. Seiler)

Anti decubitus systems

Thermally reprocessed by RKI guidelines

- · Hyper Elastic foam core with transverse and longitudinal sections
- 2-ply
- · Ventilation channels for optimum microclimate
- Optimal pressure relief and reduction of shear forces (under 18mm/Hg print run)
- PU-Tex reference abnehmbarund washable up to 95°C + wipe disinfection
- Mattress autoclavable
- · Permeable to water vapor
- Flame-retardant according to DIN EN 597 (part 1 + 2) agreed
- . 20 to 100 kg of patient weight
- Suitable for pain therapy

Reha-Star mattress II up to stage 4 (nach Prof. Seiler)

Thermally reprocessed by RKI guidelines

- . Hyper Elastic foam core with transverse and longitudinal sections
- 3-ply
- · Ventilation channels for optimum microclimate
- Optimal pressure relief and reduction of shear forces (under 18mm/Hg print run)
- PU-Tex reference abnehmbarund washable up to 95°C + Wipe disinfection
- Mattress autoclavable
- · Permeable to water vapor
- Flame-retardant according to DIN EN 597 (part 1 + 2) agreed
- . 20 to 120 kg of patient weight
- Suitable for pain therapy



Range of colours

examples for coloured features:



blue (article-no. 90.515)



blue-purple (article-no. 90.405)



dark-purple (article-no. 90.406)



orange (article-no. 90.204)



yellow (article-no. 90.128)




natural (article-no. 90.100)



green (article-no. 90.618)

Get to know our range of colour.



Make your choice ...

- · choose modell
- choose size
- · choose equipment
- · choose the colour

Do not hesitate to contact us in case of special requests!

Padding

Thanks to individual production, we can offer you the right cushioning for each specific need.

Our skai-padding is in many colours available!

You can find a preselection on our website at: www.savi.de/skai-padding

We will gladly send you the current colour selection.



| Padding thickness skai-padding 3 mm | |
|---------------------------------------|---------------|
| Skai-padding 50 cm | 99.051 |
| Skai-padding 100 cm | 99.051.100 |
| Skai-padding with vinyl windows | |
| Motifpadding Greta, 50 cm | 99.050GRETA |
| Motifpadding Paul, 50 cm | 99.050PAUL |
| Motifpadding Pauline, 50 cm | 99.050PAULINE |
| Fixed skai-padding 50 cm | 99.051.FEST |
| Fixed skai-padding with vinyl windows | 99.051.FEST |
| Eade padding | 99.052 |

Great paddings that fit really well

Wood colours

Our beds are available in many wood colours. As an alternative to "natural oiled", the beech wood is treated with the desired colour glaze — all possibilities are biologically neutral and free of harmful substances according to DIN EN 71-3 (safety of toys)!

As wood is a natural raw material, every bed looks different. Your bed will be a real unicum!



Motifpadding

Our motifpadding is easy with cords at the bars of the bed fixed.



Padding

Bound skai-padding

Our skai-padding variants includes all listed below!

Skai-padding ...

... with Velcro-fixing, height adjustmend possible! The Velcro webs are attached to the sticks.

... is using cords on the side Sticks attached. Due to the positioning of the cords always flexible giben the optimal fixation.





... for fixing with adhesive Velcro.











Fixed skai-padding

Fixed skai-padding, padding is screwed at the cot, cot-dimensions adapted so that you can use the adjustment. Attention: can only be retrofitted at Little Mila!



Fixed skai-padding with vinyl windows

The skai-padding is framing a vinyl-window. You and your child have a perfect view.





Edge padding

The padding around the edges to protect against injury.



Fixed skai-padding with edge padding

Combine the fixed skai-padding with the edge padding.



Sticks- and Frame-skaipadding

Example pictures: please see page 35. All our children's care beds are available with special versions of skai-padding.

create your own savi®-children's bed!

coloured sticks children's name:

| Create your children's bed for rehabilitation! his is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's. ust make a cross and decide how your bed should look like: Lying area of the bed? 65 x 155 cm | 13. Running gear? wheels without wheels | |
|--|--|--|
| | 13. Running gear? | |
| | | Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent. Please use this statements to create an offer about a special model to ful I the handicaps exactly. |
| | 0. Do you need perspex? Yes. Where?: No | place, date stamp, sign |
| 1. Kind of wood/wood-colour? Wood: beech pine Vood-colour: natural basalt wenge barrique teak 2. Colour feature? blue-purple dark-purple blue | Just send us your order via e Mail or use our online-form at: www.savi.de SAVI GmbH Reha Solutions Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70 | |
| green white stained white lacquered other colour: | D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60 | |

Your dealer:



SAVI GmbH Reha Solutions

Lise-Meitner-Straße 10

D-48691 Vreden

Fon: +49 (0) 25 64 - 93 09 - 70

Fax: +49(0)2564-9309-60

www.savi.de info@savi.de

Facebook: SAVI.GmbH