

Special beds for unique children – our Tom can be designed variably according to the needs of your child and therefore always offers the right solution.

Thanks to the four-winged folding sliding door, the barred sides and the flexible height adjustment possibilities, our Tom provides an optimal therapeutic support. Due to the low entry height of a little about 20 cm even smaller children can get into the bed on their own.



low entry height of about 20 cm possible,







EN 60601-2-52 DIN EN 32623 According to Medical Devices Act CE-Markering Declaration of Conformity

Best security for your child, thanks to optimal locking systems

Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- Comfortable, back-friendly height electric
- · Many bed sizes
- · Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- · Colourful padding for maximum protection
- Short delivery times

Tom' individual strengths


- . Four-winged door to open the complete front
- Castors
- · High barred sides (option)
- . Low entry height of 20 cm minimum (optional)
- Trendelenburg-/Antitrendelenburg-function (option)
- · Fixed skai-padding (also with vinyl windows) (option)



Choice of colours (see on page 68)


A lot of safety, even for particularly active children







Barred side height 135 cm, e.g. in blue-purple colour accent with wooden castors



Technical data

	The Mark Williams Indiana		- Constitution of the Cons	
Lying area	External dimension		Article-no.	
100 x 200 cm	115,5 x 223 cm		155,34,008	
90 x 200 cm	105,5 x 223 cm		155.34.008.L	
90 x 170 cm	105,5 x 193 cm		155.34.008.M	
65 x 155 cm	80,5 x 178 cm		155.34.008.5	
70 x 140 cm	85,5 x 163 cm		155.34.008.XS	
Doors				
Four-winged do	or to open the complete	front		
Four-winged do	or on both sides		99.750	+
270° door, doo	r outwardly opening		99.704	+
270° door, doo	r outwardly opening on b	ooth sides	99.704BEID	+
Locking mecha	nism			
Child-proof SAV	/Flock I			•
Child-proof SAV	/Flock I metal		40100003	+
SAVI-lock II			99.033	+
SAVI-lock III			99.034	+
SAVI-lock IV			99.034A	+
Lying height (to	p of slatted bed base)			
	ectric height adjustment			
Barred side hei		Overall heigh	i	7.17.70
72 cm		114 cm		+
113 cm		155 cm		
135 cm		177 cm		+
Special sizes on	request	177 311		+
Slatted bed bas				
No adjustment				+
Manual adjustr	Jane 1961			
Electric adjustment			99.200	+
Trendelenburg-/Antitrendelenburg-function (sloping positi		tion (sloping position)	99.211	4
Running gear		mand and sure I received a	70000	
	ning gear with total lock	ing mechanism		
Designer wood			99,885	n,a.c.
Without castors and low entrance height of 20 cm		99.883	+	
	n low entrance height of		99.881	+
Kind of wood	and an arrangement		23,000	
C. COLORS SERVICE SERV	atural oiled" (further colo	ours see p. 69)		
Mattresses	William Control	W 1/4		
Polyfoam mattr	ess			
Cold foam mattress		99.031.5	+	
AND DESCRIPTION OF THE PARTY OF	ess plus (up to stage 2)		99.031.3.PLUS	+
	ress II (up to stage 4)		99.031.3.ZWEI	+
Incontinence co			99.031.3ZA50	+
		000000000000000000000000000000000000000		
= standard	o.r. = on request	+ = accessories	n.a.c. = no additional charge	

Technical data

Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	n.a.c.
Orange	90,204	n.a.c.
Blue-purple	90,405	
Dark-purple	90.406	n.a.c.
Blue	90.515	n.a.c.
Green	90.618	n.a.c.
Equipment		
Children's name on side of bed		n.a.c.
Coloured sticks		n.a.c.
Perspexs instead of sticks		•
Sticks instead of perspexs	99.815	n.a.c.
Touch protection for height adjustment	99.900	+
External lock on side of bed	99.706	+
Feed tube hole	99.840	Q.F.
Raising aid	99.802	+
Net-roof	99.800	+
Adaption of the bed for sloping ceiling		7.0
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		
Skai-padding 50 cm height	99.051	+
Skai-padding 100 on height	99.051.100	+
Fixed skai-padding	99.051.FEST	O.F.
Fixed skai-padding with vinyl windows	99.051.FEST	0.1.
Egde padding	99.052	.+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Our skai-padding is in many colours available! You can find a pre-selection on our website at: www.savi.de/skai-padding



Thanks to its folding-sliding mechanism, the door offers the maximum stability and precise operation in every position. The design is a harmonious combination of quality and function. The entire front side can be opened in this way.

Tom adapts perfectly to the situation and makes everyday life with your child easier.

Naturally by request we can also offer customised additional doors, and custom-made sloping roofs are also not a problem for us.

Safety first – SAVI door locks

- Child-proof SAVI-lock I (base)
- Child-proof SAVI-lock I made of metal
- 3 SAVI-lock II
- SAVI-lock III: The labyrinth guide prevents the child opening the bed by accident. This lock is keyless and has no magnets. Can also be retrofitted.
- SAVI-lock IV: Here too the bed is opened by means of a labyrinth guide, but in addition this lock also operates against resistance. This lock is keyless, without magnets and can be retrofitted.
- External lock on side of bed (increases the stability and provides more secure)

Great solutions for door locks - highest safety first



Little tip:

In this catalogue you can discover many bed variants that we have them by the "ID" below the


Key data of Tom

- · Barred side height 72 cm, Barred side height 113 cm, Barred side height 135 cm / special sizes on request
- . Slatted bed base 34 104 cm or Slatted bed base 20 - 90 cm (without wheels)
- Trendelenburg-/Antitrendelenburg-function (option)
- · Choise between perspexs or sticks

270° door to fold back completely

- · Easy to use manually
- . Doors on the head and foot elements fold back completely
- 270° version can also be supplied with firm skai-padding





All functions of the bed can always be used without removing the padding beforehand!

Optimum protection thanks to fixed Skai-padding

Your child is very active? With the firmly integrated skai-padding (75 mm height), you offer your child optimum protection.

For more information on the different padding options, see page 69/70!

Optimal view thanks to integrated vinyl windows

The integrated vinyl windows fit flush with the firm skai-padding and thus offer a similar padding effect to the permanently integrated foam padding. The windows offer an optimal view. So you keep an eye on your child easily and your child feels visibly more comfortable.







Examples










ID: 20048









Tom Paddy

Tom Paddy's individual strengths

- · Low entry height of 20 cm minimum
- High barred sides 113 cm
- · Fixed skai-padding with vinyl windows
- · Castors (option)
- Trendelenburg-/Antitrendelenburg-function (option)







Now entry height of about 20 cm possible

Tom Paddy


Technical data

Lying area	External dimension	Article-no.	
90 x 200 cm 120 x 238 cm		155.34.008.L PADDY	
90 x 170 cm 120 x 208 cm		155.34.008.M PADDY	
Doors			
Four-winged do	or to open the complete front		
Four-winged do	or on both sides	99.750	+
270° door, door outwardly opening		99.704	+
270° door, door outwardly opening on both sides		99.704BEID	+
Locking mechan	nism		
Child-proof SA\			•
Child-proof SAN	/Flock I metal	40100003	+
SAVI-lock II		99.033	+
SAVI-lock III		99.034	+
SAVI-lock IV		99.034A	+
Lying height (to	p of slatted bed base)		
20 - 90 am elec	tric height adjustment		
Door height	Overall height		
123 cm	149 cm		+
Slatted bed bas	e		
Manual adjustn	nent		
Electric adjustment		99.200	+
Trendelenburg-	/Antitrendelenburg-function (sloping position)	99.211	+
Running gear			
Without castors and low entrance height of 20 cm		99.883	+
	n low entrance height of 27 cm	99.881	+
Kind of wood			
	atural oiled" (further colours see p. 69)		
Mattresses	Paris 1 20 10		
Polyfoam mattr	ess :		•
Cold foam mat		99.031.5	+
	ess plus (up to stage 2)	99.031.3.PLUS	+
Reha-star mattress II (up to stage 4)		99.031.3.ZWEI	+
Incontinence cover		99.031.3ZA50	+
= standard	o.r. = on request += accessories		

Technical data

Equipment		
External lock on side of bed	99.706	+
Raising aid	99.802	+
Net-roof	99.800	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Bluetooth-control		
Padding		
Fixed skai-padding with vinyl windows	99.051.FEST	
• = standard o.r. = on request + :	= accessories	

Fixed padding in creme or sky-blue.



create your own savie-children's bed!

Ocoloured sticks Ochildren's name:

Create your children's bed for rehabilitation!	13. Running gear? wheels without wheels		
This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's. Just make a cross and decide how your bed should look like: 1. Lying area of the bed? 65 x 155 cm	14. Do you want a mattress? Yes No 15. Do you want a mattress? Yes No 15. Do you want a mattress? Yes No 15. Do you want a mattress? Yes No 16 yes, which one? polyfoam mattress interior spring mattress cold foam mattress Reha-Star II Reha-Star Plus intercontinence cover 16. Do you need a padding? Yes No 17 yes, which one? skai-padding fixed skai-padding skai-padding for the top edge of the cot motif-padding: Paul Greta height of the padding: 17. Do you need Vinyl-windows? Yes. Where? No 18. Equipment? touch-protection net-roof raising aid motor cover (Knut) mattress-support side up aid feed tube hole removeable posts (Niklas) feed tube hole (Lexan)		
8. How many doors do you need? front side front and back sides head and footelement 9. Locking mechanism? SAVI-lock I SAVI-lock II SAVI-lock III	to the greatest possible extent. Please use this statements to create an offer about a special model to ful I the handicaps exactly.		
10. Do you need perspex? Yes. Where?: No	place, date stamp, sign		
11. Kind of wood/wood-colour? Wood: beech pine Wood-colour: natural basalt wenge barrique teak 12. Colour feature?	Very easy: Just send us your order via eMail or use our online-form at: www.savi.de SAVI GmbH Reha Solutions		
Onatural Oyellow Orange Oblue-purple Odark-purple Oblue Ogreen Owhite stained Owhite lacquered Oother colour:	Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70 D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60		