Curious how the child care bed Tom (90 x 170 cm, 113 cm barred side height) looks in your nursery?



Simply scan the QR code and place the bed visually in your nursery!



Special beds for unique children — our Tom can be designed variably according to the needs of your child and therefore always offers the right solution.

Thanks to the four-winged folding sliding door, the barred sides and the flexible height adjustment possibilities, our Tom provides an optimal therapeutic support. Due to the low entry height of a little about 20 cm even smaller children can get into the bed on their own.

low entry height of about 20 cm possible





Discover our improved corner posts in the ball chain test!



Scan the QR code and find out more about our improved corner posts (for all our beds)!



EN 60601-2-52 DIN EN 50637 According to Medical Device Regulation CE-Markering Declaration of Conformity

Best security for your child, thanks to optimal locking systems

Why a bed from SAVI?

- · Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- · Stability and reliability
- Comfortable, back-friendly height electric
- Many bed sizes
- Individual equipment options
- · Optimal eye contact with your child
- · Large, safe indoor area
- Colourful padding for maximum protection
- · Short delivery times

Tom's individual strengths


- · Four-winged folding sliding door with comfort handle
- Castors
- High barred sides (option)
- . Low entry height of about 20 cm (option)
- Trendelenburg-/Antitrendelenburg-function (option)
- Fixed Skai-padding (also with integrated vinyl windows) (option)



Choice of colours (see on page 78)

A lot of safety, even for particularly active children







Barred side height 135 cm, e.g. in blue-purple colour accent with wooden castors

Barred side height 72 cm, e.g. in blue-purple colour accent with wooden castors

Technical data

Lying area		External dime	ension	Article-no.	
100 x 200 cm		115,5 x 223	cm	155.34.008	
90 x 200 cm		105,5 x 223	cm	155.34.008.L	
90 x 170 cm		105,5 x 193	cm	155.34.008.M	
65 x 155 cm		80,5 x 178	cm	155.34.008.5	
70 x 140 cm		85,5 x 163	cm	155.34.008.XS	
Doors					1
Four-winged fol	lding sliding door				•
	lding sliding door on bo	th sides		99.750	+
	the 270° version			99.704	+
Folding door in	the 270° version on bo	th sides		99.704BEID	+
Locking mechan	nism				
Child-proof SAV	/I-lock I				•
Child-proof SAV	/I-lock I metal			40100003	+
SAVI-lock II				99.033	+
SAVI-lock III				99.034	+
SAVI-lock IV				99.034A	+
Lying height (to	p of slatted bed base)				
	ectric height adjustment	8			
Barred side heig	ght	Overall heigh	t		
72 cm		114 cm			+
113 cm		155 cm			
135 cm		177 cm			+
Special sizes on	request				+
Slatted bed bas	e				
No adjustment	(fixed)				+
Manual adjustn	nent			99.030	•
Electric adjustm	ent			99.200	+
Trendelenburg-	/Antitrendelenburg-fun	ction (tilting of the entir	e bed)	99.211	+
Running gear			ALMINA L		
Four castor runn	ning gear with total lock	king mechanism			•
Designer woode	en castors			99.885	n.a.c.
Without castors	and low entrance helg	ht of 20 cm		99.883	+
Small castors ar	nd low entrance height	of 27 cm		99.881	+
Kind of wood					
Beechwood "na	atural oiled" (further col	lours see p. 79)			•
Mattresses					
Polyfoam mattr			1	99.031.4	•
Cold foam matt	ress			99.031.5	+
Reha-star mattr	ess plus (up to stage 2)			99.031.3.PLUS	+
Reha-star mattr	ress II (up to stage 4)			99.031.3.ZWEI	+
Incontinence co	wer			99.031.3ZA50	+
• = standard	o.r. = on request	+ = accessories	n.a.c. = no addi	tional charge	

Technical data

and the second s		The second second
Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	n.a.c.
Orange	90,204	n.a.c.
Blue-purple	90.405	
Purple	90.406	n.a.c.
Blue	90.515	n.a.c.
Green (for more colours see page 78)	90.618	n.a.c.
Equipment .		
Children's name on side of bed	Kindemame	n.a.c.
Coloured sticks		n.a.c.
Perspex Instead of sticks		•
Sticks instead of perspex	99.815	n.a.c.
Touch protection for height adjustment	99.900	+
External lock on side of bed	99.706	+
Tube feed-through	99.840	O.f.
Raising aid	99.802	+
Net-roof	99.800	+
Adaption of the bed for sloping ceiling		1.0
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		
loose Skai-padding, 50 cm height	99.051	+
loose Skal-padding, 100 cm height	99.051.100	+
Fixed Skai-padding	99.051.FEST	O.f.
Fixed Skal-padding with integrated vinyl window	vs 99.051.FEST	0.f.
Egde padding	99.053	+
 = standard o.r. = on request += access 	sories n.a.c. = no additio	onal charge

Our Skai-padding is in many colours available! You can find a pre-selection on our website at: www.savi.de/en/bed-padding/



Open side front thanks to four-winged folding sliding door

Thanks to its folding sliding mechanism, the door offers the maximum stability and precise operation in every position. The design is a harmonious combination of quality and function. The entire front side can be opened in this way.

Tom adapts perfectly to the situation and makes everyday life with your child easier.

Naturally by request we can also offer customised additional doors, and custom-made sloping roofs are also not a problem for us.

Safety first – SAVI door locks

- Child-proof SAVI-lock I (base)
- Child-proof SAVI-lock I made of metal
- 3 SAVI-lock II
- SAVI-lock III: The labyrinth guide prevents the child opening the bed by accident. This lock is keyless and has no magnets. Can also be retrofitted.
- SAVI-lock IV: Here too the bed is opened by means of a labyrinth guide, but in addition this lock also operates against resistance. This lock is keyless, without magnets and can be retrofitted.
- External lock on side of bed (increases the stability and provides more secure)

Gireat solutions for door locks - highest safety first



Key data of Tom

- Barred side height 72 cm,
 Barred side height 113 cm,
 Barred side height 135 cm,
 special sizes on request
- Lying height 34 104 cm or Lying height 20 - 90 cm (without castors)
- Trendelenburg-/Antitrendelenburg-function (option)
- · Choose between perspex or sticks

270° door to fold back completely

- · Easy to use manually
- · Doors on the head and foot elements fold back completely
- 270° version can also be supplied with fixed Skai-padding







All functions of the bed can always be used without removing the padding beforehand!

Optimum protection thanks to fixed Skai-padding

Your child is very active? With the fixed Skai-padding (75 mm thickness), you offer your child optimum protection.

For more information on the different padding options, see page 79/801

Optimal view thanks to integrated vinyl windows

The integrated viryl windows fit flush with the fixed Skai-padding and thus offers a similar padding effect to the permanently integrated foam padding.

The windows offer an optimal view. So you keep an eye on your child easily and your child feels visibly more comfortable.



ID: 36009 Supply opening (hole) with closing mechanism





Scan the QR code to find out more about our fixed Skai-padding with integrated viryl windows and its padding capability.



Examples





270° door on both sides
















Tom Paddy

Tom Paddy's individual strengths

- . Low entry height of 20 cm minimum
- · High barred sides 113 cm
- · Fixed Skai-padding with integrated vinyl windows
- · Castors (option)
- Trendelenburg-/Antitrendelenburg-function (option)







Tom Paddy

Technical data

Lying area	External dimension	Article-no.	
90 x 200 cm	120 x 238 cm	155.34.008.L PADD	Y
90 x 170 cm	120 x 208 cm	155.34.008.M PADI	ΟY
Doors			
Four-winged folding sliding	door	Specifical Ed. Sp	•
Four-winged folding sliding	99.750	+	
Folding door in the 270° ve	99.704	+	
Folding door in the 270° ve	99.704BEID	+	
Locking mechanism			
Child-proof SAVI-lock I			•
Child-proof SAVI-lock I met	al	40100003	+
SAVI-lock II		99.033	+
SAVI-lock III		99.034	+
SAVI-lock IV		99.034A	+
Lying height (top of slatted	bed base)		
20 - 90 cm electric height a	adjustment (without castors)		398.6
Door height	Overall height		
113 cm	149 cm		+
Slatted bed base			
Manual adjustment			•
Electric adjustment		99.200	+
Trendelenburg-/Antitrende	99.211	+	
Running gear			
Without castors and low en	99.883	+	
Small castors an low entran	99.881	+	
Kind of wood			
Beechwood "natural oiled"	(further colours see p. 79)		•
Mattresses	21 0-12		
Polyfoam mattress		99.031.4	•
Cold foam mattress	99.031.5	+	
Reha-star mattress plus (up	99.031.3.PLUS	+	
Reha-star mattress II (up to	99.031.3.ZWEI	+	
		00 034 37450	20.45
Incontinence cover		99.031.3ZA50	+


Technical data

Equipment		
External lock on side of bed	99.706	+
Raising aid	99.802	+
Net-roof	99.800	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		

Fixed Skal-padding with integrated vinyl windows 99.051.FEST

= standard o.r. = on request + = accessories

Fixed padding in creme or light-blue.



Create your own

coloured sticks

children's name:



Create your children's bed for rehabilitation! Castors Without castors 13. Running gear? This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes Yes No 14. Do you want drawer? into consideration and make an offer about one of our bed's. Yes No 15. Do you want a mattress? Just make a cross and decide how your bed should look like: If yes, which one? opolyfoam mattress cold foam mattress interior spring mattress 1. Lying area of the bed? Reha-Star II Reha-Star Plus incontinence cover ○ 65 x 155 cm 70 x 140 cm 90 x 170 cm 90 x 200 cm 100 x 200 cm other size: 16. Do you need a padding? Yes No If yes, which one? 2. Which barred side do you need? Height (cm): fixed Skai-padding motifpadding MALU O loose Skai-padding edge-padding 3. Do you use a patient-lifter? colour of the padding: height of the padding: 4. Which lying height (top of slatted base) do you need? 17. Do you need vinyl windows? ○ 40 cm ○ 50 cm ○ 60 cm ○ 70 cm other height: Yes, Where? 5. Do you need an electric height adjustment? 18. Equipment? 6. Which slatted bed base do you need? removable posts (Niklas) raising aid lockable handset box net-roof infusion holder mattress holder tube feed-through (gap) drawers no adjustment manual adjustment step-protection for the bed base touch-protection for height adjustment foot cut-out (e.g. for patient lift) electric adjustment Trendelenburg-function central locking device for 4-roller chassis (Knut a+b only) access aid made of artificial leather 7. What kind of doors do you want? supply opening (hole) with dosing mechanism 4-winged folding sliding door 270° door side rail "hide away" folding door O lowerable side rail O electric lowerable side rail Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to 8. How many doors do you need? the greatest possible extent. front side front and back sides head and footelement Please use this statements to create an offer about a special model to fulfil the handicaps exactly. 9. Locking mechanism? SAVI-lock III SAVI-lock IV SAVI-lock I SAVI-lock II 10. Do you need perspex? O No Yes. Where?: place, date stamp, sign 11. Kind of wood/wood-colour? Wood: beech opine) wenge barrique Wood-colour: natural () basalt O teak O walnut Very easy: Just send us your order via email white stained white lacquered or use our online-form at: www.savi.de/en/home/ SAVI GmbH Reha Solutions 12. Colour feature? Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70 Orange Ored Ogrey Oblue-purple natural vellow D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60 other colour: O blue O white