

Knut



Great colour?
It's available for each of
our beds and even
more wood colours!
See page 69.



↖ e.g. in
wood colour
„Barrique“

Cuddle safely in Knut – the bed with a high degree of safety.

Thanks to the four-winged folding sliding door, the barred sides and the flexible height and inclination adjustment possibilities, this bed provides very good therapeutic support.



Trendelenburg- /Antitrendelenburg-
function



Electric height
adjustment



Head- & footelement manual or
electric adjustment (option)



More info about
Knut:
www.savi.de/beds/knut

EN 60601-2-52
DIN EN 32623

According to Medical Devices Act
CE-Markering
Declaration of Conformity

Knut

Why a bed from SAVI?

- Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- Stability and reliability
- Comfortable, back-friendly height electric
- Many bed sizes
- Individual equipment options
- Optimal eye contact with your child
- Large, safe indoor area
- Colourful padding for maximum protection
- Short delivery times



Knut's individual strengths

- Barred side 64 - 89 cm
- Trendelenburg- /Antitrendelenburg-function
- Drawer / stowage (option)
- Doors or adjustable side rails (option)
- Low entry height of 21 cm minimum (option)



All models are available with alternative side rail heights

Knut (a)

Barred side height 64 cm
lying height 41 - 76 cm or
31 - 66 cm (small castors)



Knut (b)

Barred side height 78 cm
lying height 41 - 91 cm or
31 - 81 cm (small castors)



Knut (c)

Barred side height 89 cm
lying height 31 - 101 cm or
21 - 91 cm (small castors)



* Indication

With adjustable side rails the lying height changes to 41 - 111 cm!

Technical data

Lying area	External dimension	Article-no.
100 x 200 cm	113 x 218 cm	
90 x 200 cm (size L)	103 x 218 cm	
90 x 170 cm (size M)	103 x 188 cm	
70 x 140 cm (size XS)	83 x 158 cm (Knut (c) 103 x 158 cm)	
65 x 155 cm (size S)	78 x 173 cm (Knut (c) 98 x 173 cm)	
Doors		
4-winged door to open the complete front		•
Four-winged door on both sides	99.750	+
270° door, door outwardly opening	99.704	+
270° door, door outwardly opening on both sides	99.704BEID	+
Locking mechanism		
Child-proof SAVI-lock I		•
Child-proof SAVI-lock I metal	40100003	+
SAVI-lock II	99.033	+
SAVI-lock III	99.034	+
SAVI-lock IV	99.034A	+
Lying height		
Top of slatted bed base	Barred side height	
(a) 41- 76 cm electric	64 cm ¹²⁾	• Knut (a)
(b) 41- 91 cm electric	78 cm ¹²⁾	• Knut (b)
(c) 31-101 cm electric	89 cm ¹²⁾	• Knut (c)
Additional charge for 1 adjustable side-rail*		
(+14 cm external length) manual	64 cm	99.710 +
electric		99.710EL +
manual	78 cm	99.712 +
electric		99.712EL +
manual	89 cm	99.714 +
electric		99.714EL +
Additional charge for 2 adjustable side-rails*		
(+14 cm external length) manual	64 cm	99.711 +
electric		99.711EL +
manual	78 cm	99.713 +
electric		99.713EL +
manual	89 cm	99.715 +
electric		99.715EL +
Slatted bed base		
Manual adjustment	99.030	•
Electric adjustment	99.200	+
Trendelenburg- / Antitrendelenburg-function (sloping position)		•

• = standard o.r. = on request + = accessories

(1) Indication

In combination with adjustable side rails not possible.

(2) Indication

All models are available with alternative side rail heights. We advise you gladly.

Technical data

Running gear		
Four castor running gear with total locking mechanism		•
Designer wooden castors	99.885	n.a.c.
Without or small rolls ¹³⁾ (21-91 cm)	99.882	n.a.c.
Kind of wood		
Beechwood „natural oiled“ (further colours see p. 69)		•
Mattresses		
Polyfoam mattress		•
Cold foam mattress	99.031.5	o.r.
Reha-star mattress plus (up to stage 2)	99.031.3.PLUS	+
Reha-star mattress II (up to stage 4)	99.031.3.ZWEI	+
Incontinence cover	99.031.3ZA50	+
Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	n.a.c.
Orange	90.204	•
Blue-purple	90.405	n.a.c.
Dark-purple	90.406	n.a.c.
Blue	90.515	n.a.c.
Green	90.618	n.a.c.
Equipment		
Children's name on side of bed		n.a.c.
Coloured sticks		n.a.c.
Perspers instead of sticks		o.r.
Perspers in the lower area of the head/foot end (standard on C)		o.r.
External lock on side of bed	99.706	+
Feed tube hole	99.840	o.r.
Raising aid	99.802	+
Sit up aid	99.820	+
2 drawers	99.040	+
Motor cover	99.803	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		
Motifpadding Greta, 50 cm	99.050GRETA	+
Motifpadding Paul, 50 cm	99.050PAUL	+
Motifpadding Pauline, 50 cm	99.050PAULINE	+
Skai-padding 50 cm height	99.051	+
Fixed skai-padding	99.051.FEST	o.r.
Fixed skai-padding with vinyl windows	99.051.FEST	o.r.
Edge padding	99.052	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

* Indication

If you choose KNUT C with drop down sides (artikel-no.: 99.710 - 99.715) the lying height will change to 41 - 111 cm.

Open side front thanks to four-winged sliding door

Thanks to its folding-sliding mechanism, the door offers the maximum stability and precise operation in every position. The entire front side can be opened in this way. The design is a harmonious combination of quality and function. Naturally by request we can also offer customised additional doors, and custom-made sloping roofs are also not a problem for us.



ID: 19118



ID: 20218



ID: 20163

Key data of Knut

- Barred side height 64 cm,
Barred side height 78 cm,
Barred side height 89 cm,
in all model variants or different
height, on request
- Lying height

Knut (a)

41 - 76 cm
(with small castors 31 - 66 cm)

Knut (b)

41 bis 91 cm
(with small castors 31 - 81 cm)

Knut (c)

31 bis 101 cm
(with small castors 21 - 91 cm)

Hinweis: With adjustable side rails the lying height changes to 41 cm!


- Trendelenburg-/Antitrendelenburg-function
- Perspex (option)
- Four-winged sliding door
- Adjustable side rail (option)



ID: 18852



ID: 19842



ID: 19463

✓ The plus when it comes to intensive care

Tubing furnishing

Our tubing furnishing is available as follows:



tubing with sliding perspex (recloseable)



tubing on the middle foot element



ID: 17139



ID: 14765



ID: 17018



ID: 20315



ID: 17496

Adjustable side rails



ID: 20355



ID: 14549



ID: 17451

Recessed safety grip

- Easy to use manually
- The safety gap on both recesses prevents unsupervised opening
- Minimal operating strength required (slider buttons)
- Grid adjustment



ID: 20350





Versions adjustable side rails

Our adjustable side rail is available as follows:



ID: 18688

Version 1

Adjustable side rail with 2 top elements (choise manually or electrically adjustable).



ID: 19926

Version 2

Adjustable side rail with 1 top element. You can remove the top element if you drop down the side rail.



ID: 18568

Special version 3

Adjustable side rail with hinged top element.

Attention: This version is only possible in combination with an electric side rail!

Benefits of electric adjustable sides

- Easy to use at the push of a button
- Lockable control element prevents unsupervised opening
- Can also be operated by people with disabilities
- Minimal operating strength required (pushbutton)
- Infinitely adjustable



ID: 19046



Knut

270° door to fold back completely

- Easy to use manually
- Doors on the head and foot elements fold back completely
- 270° version can also be supplied with firm skai-padding



„These are beastly good beds!“



ID: 16796



ID: 16017

Examples

Indication

Our children's therapy beds come in a choice of white painted or white stained (wood grain visible) finishes.



ID: 13623



ID: 19333



ID: 18732



ID: 19489

Knut

Optimum protection thanks to fixed skai-padding

Your child is very active? With the firmly integrated skai-padding (75 mm height), you offer your child optimum protection.

For more information on the different padding options, see page 69/ 70!

Our skai-padding is in many colours available!

You can find a pre-selection on our website at:

www.savi.de/skai-padding



ID: 20506



ID: 18581



ID: 16277



ID: 16923



ID: 18323



Optimal view thanks to integrated vinyl windows

The integrated vinyl windows fit flush with the firm skai-padding and thus offer a similar padding effect to the permanently integrated foam padding. The windows offer an optimal view. So you keep an eye on your child easily and your child feels visibly more comfortable.



Tubing furnishing

On request, we manufacture all care beds, with tubing furnishing.




Edge padding



Stick- / Framepadding



Knut (b) 40



Knut (b) 40's individual strengths

- 64 cm barred side height, by removable push-on element (40 cm barred height)
- Trendelenburg-/Antitrendelenburg-function
- Drawer/stowage (option)
- Adjustable barred sides on both sides (option)

Barred side heights

- barred side height 40 cm + push-on element 24 cm (total height 64 cm to comply with legal requirements)
- Slatted bed base 41 - 91 cm

variable and safe



adjustable side rails
on both sides (option)



low side rails
due to removable
push-on element
on both sides (option)



one side adjustable barred side
with 24 cm
push-on element
(standard)

Technical data

Lying area		External dimension	
90 x 200 cm (size L)		106,5 x 232 cm	
90 x 170 cm (size M)		106,5 x 202 cm	
Lying height		Barred side height	
Top of slatted bed base		64 cm bzw. 40 cm	
41 - 91 cm electric		64 cm bzw. 40 cm	
Adjustable side rail		Bar height	
(+14 cm external length)	manual	64 cm	99.710 •
	electric	64 cm	99.710EL +
2 adjustable side rails		Bar height	
(+14 cm external length)	manual	64 cm	99.711 +
	electric	64 cm	99.711EL +
Slatted bed base			
Manual adjustment		99.030 •	
Electric adjustment		99.200 +	
Trendelenburg- / Antitrendelenburg-function (sloping position)		•	

• = standard o.r. = on request + = accessories



Technical data

Kind of wood		
Beechwood „natural oiled“ (further colours see p. 69)		•
Mattresses		
Polyfoam mattress		•
Cold foam mattress	99.031.5	+
Reha-star mattress plus (up to stage 2)	99.031.3.PLUS	+
Reha-star mattress II (up to stage 4)	99.031.3.ZWEI	+
Incontinence cover	99.031.3ZA50	+
Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	n.a.c.
Orange	90.204	•
Blue-purple	90.405	n.a.c.
Dark-purple	90.406	n.a.c.
Blue	90.515	n.a.c.
Green	90.618	n.a.c.
Equipment		
Children's name on side of bed		n.a.c.
Coloured sticks		n.a.c.
Perspex instead of sticks		o.r.
Perspex in the lower area of the head / foot end (standard on C)		o.r.
Feed tube hole	99.840	o.r.
Raising aid	99.802	+
Sit up aid	99.820	+
2 drawers	99.040	+
Motor cover	99.803	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Bluetooth control		
Padding		
Motifpadding Greta, 50 cm	99.050GRETA	+
Motifpadding Paul, 50 cm	99.050PAUL	+
Motifpadding Pauline, 50 cm	99.050PAULINE	+
Skai-padding 50 cm height	99.051	+
Fixed skai-padding	99.051.FEST	o.r.
Fixed skai-padding with vinyl windows	99.051.FEST	o.r.
Edge padding	99.052	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Create your own SAVI®-children's bed!

Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our bed's.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

- ☐ 65 x 155 cm ☐ 70 x 140 cm ☐ 90 x 170 cm
☐ 90 x 200 cm ☐ 100 x 200 cm other size: _____

2. Which barred side do you need? Height (cm): _____

3. Do you use a patient-lifter? ☐ Yes ☐ No

4. Which lying height (top of slatted base) do you need?

- ☐ 40 cm ☐ 50cm ☐ 60 cm ☐ 70 cm other height: _____

5. Do you need an electric height adjustment? ☐ Yes ☐ No

6. Which slatted bed base do you need?

- ☐ no adjustment ☐ manual adjustment ☐ protection for the bed base
☐ electric adjustment ☐ Trendelenburg function

7. What kind of doors do you want?

- ☐ 4-winged door ☐ 270° door ☐ folding door
☐ adjustable side rail ☐ electric adjustable side rail

8. How many doors do you need?

- ☐ front side ☐ front and back sides ☐ head and footelement

9. Locking mechanism?

- ☐ SAVI-lock I ☐ SAVI-lock II ☐ SAVI-lock III ☐ SAVI-lock IV

10. Do you need perspex?

- ☐ Yes. Where?: _____ ☐ No

11. Kind of wood/wood-colour? Wood:

- ☐ beech ☐ pine
 Wood-colour: ☐ natural ☐ basalt ☐ wenge ☐ barrique ☐ teak

12. Colour feature?

- ☐ natural ☐ yellow ☐ orange ☐ blue-purple ☐ dark-purple ☐ blue
☐ green ☐ white stained ☐ white lacquered ☐ other colour: _____
☐ coloured sticks ☐ children's name: _____

13. Running gear? ☐ wheels ☐ without wheels

14. Do you want containers? ☐ Yes ☐ No

15. Do you want a mattress? ☐ Yes ☐ No

If yes, which one?

- ☐ polyfoam mattress ☐ interior spring mattress ☐ cold foam mattress
☐ Reha-Star II ☐ Reha-Star Plus ☐ intercontinence cover

16. Do you need a padding? ☐ Yes ☐ No

If yes, which one?

- ☐ skai-padding ☐ fixed skai-padding ☐ skai-padding for the top edge of the cot
 motif-padding: ☐ Pauline ☐ Paul ☐ Greta

height of the padding: _____ colour of the padding: _____

17. Do you need Vinyl-windows?

- ☐ Yes. Where? _____ ☐ No

18. Equipment?

- ☐ touch-protection ☐ net-roof ☐ raising aid ☐ motor cover (Knut)
☐ mattress-support ☐ side up aid ☐ feed tube hole ☐ removeable posts (Niklas)
☐ feed tube hole (Lexan)

☐ Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

☐ Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, date

stamp, sign

Very easy:

Just send us your order via eMail
 or use our online-form at: www.savi.de

SAVI GmbH Reha Solutions

Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70
 D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60

