

Lola



Small heroes very big – Lola is the ideal bed for toddlers.

This bed is ideal for small children because it has a fixed, high basic lying height, making nursing more comfortable.

Its also easy on the back, thereby seeing to restful sleep.



More info about
Lola:
www.savi.de/beds/lola

EN 60601-2-52
DIN EN 32623

According to Medical Devices Act
CE-Markering
Declaration of Conformity

Why a bed from SAVI?

- Great safety for your child
- Biologically neutral, pollutant-free (DIN EN 71-3)
- Stability and reliability
- Comfortable, back-friendly height electric
- Many bed sizes
- Individual equipment options
- Optimal eye contact with your child
- Large, safe indoor area
- Colourful padding for maximum protection
- Short delivery times

Lola's individual strengths

- Good working height / slatted bed base 85 cm
- Cupboards (option)



Choice of colours (see on page 68)



door open
→



barred side height 65 cm
z.B. im Farbakzent grün



barred side height 65 cm and cupboards
e.g. in green colour accent

Technical data

Lying area	External dimension	Article-no.	
100 x 200 cm	113 x 206 cm		
90 x 200 cm	103 x 206 cm		
Doors			
4-winged door to open the complete front			•
Four-winged door on both sides		99.750	+
270° door, door outwardly opening		99.704	+
270° door, door outwardly opening on both sides		99.704BEID	+
Locking mechanism			
Child-proof SAVI-lock I			•
Child-proof SAVI-lock I metal		40100003	+
SAVI-lock II		99.033	+
SAVI-lock III		99.034	+
SAVI-lock IV		99.034A	+
Lying height (top of slatted bed base)			
85 cm			•
Barred side height		Overall height	
65 cm	150 cm	150.85.002	•
90 cm	175 cm	175.85.002	+
Special sizes on request			+
Slatted bed base			
Manual adjustment		99.030	+
Electric adjustment		99.200	+
Running gear			
Four castor running gear with total locking mechanism			o.r.
Designer wooden castors			o.r.
Kind of wood			
Pinewood			n.a.c.
Beechwood „natural oiled“ (further colours see p. 69)			•
Mattresses			
Polyfoam mattress		99.031.4	•
Cold foam mattress		99.031.5	+
Reha-Star Plus (Stadium 2)		99.031.3 PLUS	+
Reha-Star II (Stadium4)		99.031.3.ZWEI	+
Incontinence cover		99.031.3.ZA50	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge

Technical data

Coloured features		
Natural	90.100	n.a.c.
Yellow	90.128	•
Orange	90.204	n.a.c.
Blue-purple	90.405	n.a.c.
Dark-purple	90.406	n.a.c.
Blue	90.515	n.a.c.
Green	90.618	n.a.c.
Equipment		
Children's name on side of bed		n.a.c.
Coloured sticks to go to the coloured feature		n.a.c.
Perspex instead of sticks		o.r.
External lock on side of bed	99.706	+
Feed tube hole	99.840	o.r.
Raising aid	99.802	+
Sit up aid	99.820	+
Adaption of the bed for sloping ceiling		o.r.
Cupboards (2 pieces)	99.041	+
Mattress holder	99.891	+
Protection for the bed base in case of kicks	99.892	+
Wood care set	99.880	+
Padding		
Motifpadding Greta, 50 cm	99.050GRETA	+
Motifpadding Paul, 50 cm	99.050PAUL	+
Motifpadding Pauline, 50 cm	99.050PAULINE	+
Skai-padding 50 cm height	99.051	+
Fixed skai-padding	99.051.FEST	o.r.
Fixed skai-padding with vinyl windows	99.051.FEST	o.r.
Edge padding	99.052	+

• = standard o.r. = on request + = accessories n.a.c. = no additional charge



Lola

Key data of Lola

- barred side height of 65 - 90 cm, other heights on request
- Fixed lying height 85 cm
- Child-proof SAVI-lock
- Cupboards (option)
- Perspex (option)



ID: 17042
fixed skai-padding in white



ID: 15875
barred sides adjustable



ID: 16211
lying area 100 x 100 cm,
barred side height 40 cm



ID: 19014
head element an barred side
to open with hinged element



ID: 18804
white painted

"I like to dream
about my friends
and the beach."



ID: 19540
tubing with sliding perspex



ID: 13874
cupboards



ID: 17769
all sticks

Create your own SAVI®-children's bed!

Create your children's bed for rehabilitation!

This is a list which you can use to create your children's bed for rehabilitations. We will take your wishes into consideration and make an offer about one of our beds.

Just make a cross and decide how your bed should look like:

1. Lying area of the bed?

- 65 x 155 cm 70 x 140 cm 90 x 170 cm
 90 x 200 cm 100 x 200 cm other size: _____

2. Which barred side do you need? Height (cm): _____

3. Do you use a patient-lifter? Yes No

4. Which lying height (top of slatted base) do you need?

- 40 cm 50cm 60 cm 70 cm other height: _____

5. Do you need an electric height adjustment? Yes No

6. Which slatted bed base do you need?

- no adjustment manual adjustment protection for the bed base
 electric adjustment Trendelenburg function

7. What kind of doors do you want?

- 4-winged door 270° door folding door
 adjustable side rail electric adjustable side rail

8. How many doors do you need?

- front side front and back sides head and footelement

9. Locking mechanism?

- SAVI-lock I SAVI-lock II SAVI-lock III SAVI-lock IV

10. Do you need perspex?

- Yes. Where?: _____ No

11. Kind of wood/wood-colour? Wood: beech pine

- Wood-colour: natural basalt wenge barrique teak

12. Colour feature?

- natural yellow orange blue-purple dark-purple blue
 green white stained white lacquered other colour: _____
 coloured sticks children's name: _____

13. Running gear? wheels without wheels

14. Do you want containers? Yes No

15. Do you want a mattress? Yes No

If yes, which one?

- polyfoam mattress interior spring mattress cold foam mattress
 Reha-Star II Reha-Star Plus intercontinence cover

16. Do you need a padding? Yes No

If yes, which one?

- skai-padding fixed skai-padding skai-padding for the top edge of the cot
 motif-padding: Pauline Paul Greta

height of the padding: _____ colour of the padding: _____

17. Do you need Vinyl-windows?

- Yes. Where? _____ No

18. Equipment?

- touch-protection net-roof raising aid motor cover (Knut)
 mattress-support side up aid feed tube hole removeable posts (Niklas)
 feed tube hole (Lexan)

Please use this statements to create an offer, where you choose a model, which fulfils the handicaps to the greatest possible extent.

Please use this statements to create an offer about a special model to fulfil the handicaps exactly.

place, date

stamp, sign

Very easy:

Just send us your order via eMail
or use our online-form at: www.savi.de

SAVI GmbH Reha Solutions

Lise-Meitner-Str. 10 Fon: +49 (0) 25 64 - 93 09-70

D-48691 Vreden Fax: +49 (0) 25 64 - 93 09-60

