

adaption form



SAVI GmbH Reha Solutions
Lise-Meitner-Straße 10
D-48691 Vreden
phone: +49 (0) 25 64 93 09 70
fax: +49 (0) 25 64 93 09 60
info@savi.de

☐ offer ☐ order

company:

phone:

fax:

e-mail:

contact person:

commission:

typ of bed	doors	frame	colour	lying area	barred side
☐ Tom	☐ 4 doors	☐ slatted bed	☐ beech natural	☐ 65 x 155 cm	☐ 60 cm
☐ Tom-Paddy	☐ 8 doors	☐ base - no adjustment	☐ teak	☐ 70 x 140 cm	☐ 64 cm
☐ Knut a	☐ 270° door	☐ electrically adjustable slat frame	☐ barrique	☐ 90 x 170 cm	☐ 72 cm
☐ Knut b	☐ folding door	☐ trendelenburg-function	☐ wenge	☐ 90 x 200 cm	☐ 78 cm
☐ Knut b 40	☐ folding sliding door	☐ kick protection	☐ basalt	☐ 100 x 200 cm	☐ 89 cm
☐ Knut c	☐ 1 manually lowerable side rail (not suitable for: Savoir-Vivre, Jonas & Tom)	☐ fixed bed base	☐ white lacquered	☐ 120 x 200 cm	☐ 100 cm
☐ Niklas	☐ 2 manually lowerable side rails (not suitable for: Savoir-Vivre, Jonas & Tom)		☐ white stained	☐ 140 x 200 cm	☐ 113 cm
☐ Savoir-Vivre	☐ 1 electr. lowerable side rail (not suitable for: Savoir-Vivre, Jonas & Tom)		☐ walnut	other:	☐ 125 cm
☐ Little Mila	☐ 2 electr. lowerable side rails (not suitable for: Savoir-Vivre, Jonas & Tom)		colour accent of choice:		☐ 135 cm
☐ Lola	☐ 1 side rail "hide away" (only for Knut a+b) ☐ 40 cm ☐ 64 cm		☐ coloured sticks at:		☐ 150 cm
☐ Jonas	☐ 2 side rails "hide away" (only for Knut a+b) ☐ 40 cm ☐ 64 cm				☐ 170 cm
					other:

mattress	lying height	padding (see back)	locking/lexan (see back)	equipment
☐ polyfoam	☐ 20 cm (Tom)	☐ loose skai-padding	☐ SAVI-lock I	☐ touch protection
☐ interior spring	☐ 21 cm (Knut c)	☐ fixed skai-padding (please fill out the back)	☐ SAVI-lock II	☐ net roof
☐ cold foam	☐ 27 cm (Tom)	☐ edge padding (only for fixed padding)	☐ SAVI-lock III	☐ infusion holder
☐ Reha Star II	☐ 30 cm	☐ motif-padding	☐ SAVI-lock IV	☐ supply port
☐ Reha-Star Plus	☐ 40 cm	☐ vinyl-window (only for fixed padding)	plexiglass/lexan height:	☐ hose passage
☐ incontinence cover	☐ 50 cm	height & colour:	(For plexiglas/lexan, please fill out the back!)	☐ lifting pole
	☐ 60 cm			☐ mattress holder
	☐ 70 cm			☐ drawers
	other:			☐ castor
		(For padding and vinyl-windows, please fill out the back!)		☐ central locking device for 4-castor running gear
				☐ motor cover
				☐ additional external door lock at the bottom of the bed
				☐ entry aid
				☐ NEW: design elements on the bedside
				☐ foot cut-out (7 cm height)
				☐ patient lift cut-out (15 cm height)

Padding (please fill in the diagram below)

Fixed padding / Vinyl-window / Loose padding

Enter the appropriate letter **F**, **V** or **L** in the fields.

material:

☐ faux leather

□ Bisonyl

example:

fixed skai-padding +
vinyl-windows

fixed skai-padding

A net of a rectangular prism is shown. It consists of six rectangular faces arranged in a cross-like pattern. The faces are labeled as follows: 'back' (top face), 'head' (right face), 'foot' (bottom face), and 'front' (left face). The net is composed of several rectangles of different sizes, with the central horizontal row being the longest. The 'back' face is at the top, 'head' is on the right, 'foot' is at the bottom, and 'front' is on the left. The net is drawn with black outlines on a white background.

Perspex or wooden sticks (please fill in the diagram below)

Perspex / Sticks

Enter the appropriate letter **P** (perspex) or **S** (sticks) in the fields.

The diagram shows a 3D box with four labels: 'back' at the top right, 'head' on the right side, 'foot' at the bottom left, and 'front' at the bottom center. The box is composed of a grid of 16 squares arranged in 4 rows and 4 columns. The top row of squares is labeled 'back'. The rightmost column of squares is labeled 'head'. The bottom row of squares is labeled 'foot'. The bottom-left square is labeled 'front'. The remaining squares are empty.

comments:

--

date:

signature: